

sleep•safe harmony

Instructions for Use

Software version: 2.3

Edition: 7A-2015

Part no.: M56 482 1

CE 0123



**FRESENIUS
MEDICAL CARE**

Table of contents

1 Index

2 Important information

2.1	How to use the Instructions for Use	2-1
2.2	Significance of warnings	2-3
2.3	Significance of notes	2-3
2.4	Significance of tips.....	2-3
2.5	Brief description.....	2-3
2.6	Intended use	2-4
2.6.1	Intended purpose.....	2-4
2.6.2	Specification of use	2-5
2.6.3	Side effects.....	2-5
2.6.4	Contraindications.....	2-5
2.6.5	Interactions with other systems	2-6
2.6.6	Therapy restrictions	2-6
2.6.7	Target group	2-6
2.7	Considerations for working on the device.....	2-6
2.8	Expected service life	2-7
2.9	Duties of the responsible organization	2-7
2.10	Operator responsibility	2-8
2.11	Disclaimer of liability.....	2-9
2.12	Warnings	2-9
2.12.1	Warnings, hygiene.....	2-9
2.12.2	Warnings, therapy	2-10
2.12.3	Warnings, system.....	2-10
2.12.4	Warnings, electrical	2-12
2.12.5	Warnings, consumables and accessories	2-12
2.13	Addresses	2-13

3 Design

3.1	Views	3-1
3.1.1	Front view	3-1
3.1.2	Rear view.....	3-2
3.2	User interface.....	3-3
3.2.1	Screen colors.....	3-5
3.3	General procedure for entering parameters	3-6
3.3.1	Entering text	3-6

3.3.2	Entering numbers	3-7
3.4	Selecting and editing options or parameters.....	3-8
3.4.1	Selecting and editing the solution	3-11
3.5	Profiling	3-15

4 Operation

4.1	Preparing the device	4-1
4.1.1	Functional test	4-1
4.2	Preparing for treatment.....	4-4
4.2.1	Preparing the supplies	4-4
4.2.2	Confirming the prescription.....	4-8
4.2.3	Inserting the <i>sleep•safe</i> Set.....	4-9
4.2.4	Connecting the solution bags	4-12
4.2.4.1	Connected solution bags do not match the prescription	4-15
4.2.4.2	Continuing the connection of the solution bags.....	4-16
4.2.5	Connecting the drain line system.....	4-17
4.2.6	Priming.....	4-19
4.2.6.1	Checking the fill level in the patient line	4-19
4.2.7	Confirming treatment data	4-20
4.3	Starting the treatment	4-21
4.3.1	Connecting the patient.....	4-21
4.3.2	Starting the treatment	4-24
4.3.3	Terminating the initial outflow	4-25
4.3.4	Treatment overview	4-27
4.4	Terminating the treatment	4-28
4.4.1	Disconnecting the patient	4-28
4.4.2	Treatment results.....	4-32
4.4.3	Draining the line set.....	4-36
4.4.4	Removing the <i>sleep•safe</i> Set	4-37
4.5	Special operational features of the <i>sleep•safe</i> harmony	4-39
4.5.1	Personalizing the <i>sleep•safe</i> harmony for a patient	4-39
4.5.2	Preparing the <i>sleep•safe</i> harmony for a treatment without a patient card	4-45
4.5.3	Removing the <i>sleep•safe</i> Set after premature termination of a treatment.....	4-46
4.5.4	Disconnection after premature termination of a treatment.....	4-47
4.5.4.1	Switching the <i>sleep•safe</i> harmony off	4-51
4.6	Options / changing data prior to treatment.....	4-52
4.6.1	Patient options	4-52
4.6.1.1	Selecting the language	4-52
4.6.1.2	Editing personal data.....	4-53
4.6.1.3	Editing patient parameters	4-54
4.6.1.4	Changing the additional outflow	4-57
4.6.1.5	Access levels.....	4-59
4.6.2	Therapy options	4-61
4.6.2.1	Creating a prescription	4-61
4.6.2.2	Selecting a prescription	4-64
4.6.2.3	Editing a prescription	4-65
4.6.2.4	Deleting a prescription.....	4-67
4.6.2.5	Displaying treatment reports	4-68
4.6.2.6	Editing the prescription schedule	4-70

4.6.3	Device options	4-72
4.6.3.1	Setting the date and time	4-72
4.6.3.2	LAN setup data	4-73
4.6.3.3	Relocating the <i>sleep•safe harmony</i> within a building	4-74
4.6.3.4	Setting the volume	4-79
4.6.3.5	Setting the brightness	4-80
4.6.3.6	Setting the screen saver	4-82
4.6.3.7	Displaying device information.....	4-83
4.6.3.8	Clinic login	4-85
4.6.3.9	Service login.....	4-87
4.6.3.10	Clinic / Service logout.....	4-89
4.7	Therapy options during treatment	4-91
4.7.1	Pause treatment	4-92
4.7.2	Pause with disconnection	4-94
4.7.3	Manual outflow	4-97
4.7.4	Skip phase.....	4-101
4.7.5	Terminating the treatment	4-102
4.7.6	Displaying the treatment report	4-104
4.8	Performing a PD-Plus treatment	4-104
4.8.1	Connecting the patient.....	4-104
4.8.2	Starting the PD-Plus cycle.....	4-105
4.8.2.1	Stopping the initial outflow	4-105
4.8.3	Connecting the patient for nighttime treatment.....	4-108

5 Alarm processing

5.1	Screen layout	5-1
5.2	Resetting the audible alarm.....	5-2
5.3	Color identification of the screen messages	5-3
5.4	Screen messages	5-6
5.5	Premature termination of a treatment after a system error / device fault	5-8
5.6	Emergency shutdown	5-8
5.7	Power failure	5-10
5.8	Screen failure	5-11

6 Cleaning / disinfection

6.1	Cleaning	6-1
-----	----------------	-----

7 Functional description

7.1	Description of functional procedures.....	7-1
7.2	Therapy types	7-2
7.2.1	Standard prescription	7-2
7.2.2	Adapted APD prescription	7-3

7.2.3	Tidal prescription	7-4
7.2.4	Basic prescription	7-6
7.2.5	PD-Plus prescription	7-7
7.3	Therapy options	7-8
7.3.1	Volume optimization	7-9
7.3.1.1	Permitted patient volume	7-9
7.3.1.2	Permitted residual volume	7-9
7.3.1.3	Permitted reduction of the inflow volume	7-10
7.3.2	Time optimization	7-10
7.3.2.1	Permitted reduction of the dwell duration	7-10
7.3.2.2	Catheter performance	7-11

8 Consumables, accessories, additional equipment

8.1	Consumables	8-2
8.1.1	Dialysis solutions	8-2
8.1.2	Single use items	8-3
8.1.3	Surface disinfection / surface cleaning	8-4
8.2	Accessories	8-4
8.3	Additional equipment	8-4
8.4	<i>sleep•safe harmony</i>	8-4

9 Installation

9.1	Connection requirements	9-1
9.1.1	Environment	9-1
9.1.2	Power supply (electrical power network)	9-1
9.1.2.1	General requirements	9-1
9.2	Installation requirements	9-2
9.3	Setting up the <i>sleep•safe harmony</i>	9-2
9.4	Installation after shipment / transport outside of buildings	9-3
9.4.1	Visual check after transport	9-3
9.4.2	System check	9-4
9.4.3	Extended system test	9-4

10 Transport / storage

10.1	Transport within buildings	10-1
10.2	Shipment / transport outside of buildings	10-1
10.3	Storage	10-1
10.3.1	Storage conditions	10-1
10.4	Environmental compatibility / disposal	10-2

11 Technical safety checks/maintenance procedures

11.1	Important information on the technical safety checks/maintenance procedures	11-1
------	---	------

12 Specifications

12.1	Dimensions and weight	12-1
12.2	Identification label (<i>sleep-safe harmony</i> product marking)	12-1
12.3	Electrical safety	12-2
12.4	Electrical power supply	12-2
12.5	Fuses	12-3
12.6	Information on electromagnetic compatibility	12-4
12.6.1	EMC guidance and manufacturer's declaration	12-6
12.7	Operating conditions	12-10
12.8	Storage conditions	12-10
12.9	External connection options	12-10
12.10	Batteries	12-11
12.11	Parameters	12-12
12.12	Materials used	12-14
12.13	Factory settings	12-16

13 Definitions

13.1	Definitions and terms	13-1
13.2	Abbreviations	13-1
13.3	Symbols	13-2
13.4	Certificates	13-3

14 Options

15 Appendix

15.1	<i>sleep-safe harmony</i> icons list	15-1
15.2	Disconnection with the PIN-Reload device	15-3
15.3	Connection with the PIN-Reload device	15-5
15.4	Instructions on the use of "Free software"	15-8

1 Index

A

Abbreviations 13-1
Access levels 4-59
Accessories 2-12, 8-1, 8-4
Adapted APD prescription 7-3
Additional equipment 8-4
Additional outflow 4-57
Addresses 2-13
Alarm output 12-11
Alarm processing 5-1
Appendix 14-1, 15-1
Atmospheric pressure 10-1, 12-10
Automated peritoneal dialysis (APD) 2-3
Auxiliary materials 12-16

B

Basic prescription 7-6
Batteries 12-11
Brief description 2-3

C

Catheter performance 7-11
CE mark 12-1
Certificates 13-3
Cleaning / disinfection 6-1
Considerations for working on the device 2-6
Consumables 2-12, 4-4, 4-5, 4-38, 4-47, 8-1, 8-2
Consumables, accessories, additional equipment 8-1
Contraindications 2-5

D

Damage 9-3
Date of manufacture 12-1
Defects 2-11, 9-3
Definitions 13-1
Design 3-1

Dimensions 12-1

Disclaimer of liability 2-9

Disposal 10-2

Draining the line set 4-36

Duties of the responsible organization 2-7

E

Electrical hazards 2-12

Electrical power supply 12-2

Electrical safety 12-2

Electromagnetic compatibility, information 12-4

Electromagnetic emissions 12-6

Electromagnetic radiation 9-2

EMC electromagnetic compatibility 12-4

EMC guidance and manufacturer's declaration 12-6

Environmental compatibility 10-2

Extended system test 9-4

Extension cable 2-12

External connection options 12-10

F

Front view 3-1

Functional description 7-1

Functional procedures, description 7-1

Functional test 4-1

Fuses 12-3

H

How to use the Instructions for Use 2-1

Hygiene 2-9

Hygiene practices 2-9, 4-4, 4-21, 4-28, 4-47, 4-96, 5-8, 5-11

Hygiene regulations 2-9, 4-4, 4-21, 4-28, 4-47, 4-96, 5-8, 5-11

I

Icons list 15-1

Identification label 12-1

Illustrations 2-2

Important information 2-1

Installation 9-1, 9-2

Installation requirements 9-2

Instruction 2-1

Intended purpose 2-4

Intended use 2-4

K

Keys and buttons 2-1

L

LAN 12-11

Level difference 4-1

Line set 2-10

Line voltage 12-2

M

Maintenance procedures 11-1

Materials used 12-14

Metals 12-15

N

Note symbol, significance 2-3

Numbered instruction 2-2

O

Operating conditions 2-8, 2-10, 9-2, 12-10

Operation 4-1

Operator responsibility 2-8

P

Parameters 12-12

Patient card 2-8, 3-2, 4-2, 4-4, 4-9, 4-45

Patient connector 2-9, 4-4, 4-21, 4-28, 4-47, 4-96, 5-8, 5-11, 8-3

Patient parameters 4-54

PatientCardplus 8-4

PD-Plus prescription 7-7

Personalizing 4-39

Plastics 12-14

Power cable 2-12

Power requirements 12-1

Power strip 2-12
Power supply connection 4-1, 12-2
Preparation 4-1, 4-4
Protection class 12-2

R

Relative humidity 10-1, 12-10
Reproduction 2-2
Risk of infection 2-9
Risk of injury 2-10, 2-11, 2-12

S

Screen 9-7
Serial number 12-1
Service life 2-7
Service Manual 11-1
Service support 2-11, 4-10, 9-3, 9-7, 9-9
Service support, international 2-13
Service support, local 2-13
Setting the brightness 4-80
Setting the volume 4-79
Side effects 2-5
sleep•safe Set 2-10, 2-12, 4-4
sleep•safe Set Plus 2-12
Specification of use 2-5
Specifications 12-1
Standard prescription 7-2
Status indicator 3-1
Storage 10-1
Storage conditions 10-1, 12-10
Symbols 13-2
System check 9-4

T

Target group 2-6
Technical Safety Checks 11-1
Temperature 10-1, 12-10
Temperature variations 9-4
Terms 13-1

Therapy 2-10
Therapy options 7-8
Therapy types 7-2
Tidal prescription 7-4
Time optimization 7-10
Tip symbol, significance 2-3
Transport 10-1
Transport / storage 10-1
Transport case 9-3
Transport outside of buildings 9-3, 10-1
Transport within buildings 10-1
Treatment data 4-2
Treatment duration 2-5, 4-32, 4-34, 4-69
Type identification 12-1

V

Visual check 9-3

W

Warning symbol, significance 2-3
Warnings 2-9
Warnings, consumables and accessories 2-12
Warnings, electrical 2-12
Warnings, hygiene 2-9
Warnings, system 2-10
Warnings, therapy 2-10
Weight 12-1

Y

Year of manufacture 12-1

2 Important information

2.1 How to use the Instructions for Use

Identification	The document can be identified by the following information on the title page and on the labels, if any: <ul style="list-style-type: none"> – Device software version – Technical document edition – Technical document part number 						
Footer	The footer contains the following information: <ul style="list-style-type: none"> – Company name – Device type – The English abbreviation for the document type and the international abbreviation for the document language, e.g., IFU-EN means Instructions for Use in English. – The edition, for example, 3A-2014 refers to edition 3A released in 2014 – The page identification, e.g., 1-3 refers to chapter 1, page 3 						
Organization of the chapters	To facilitate the use of documents from Fresenius Medical Care, the organization of the chapters has been standardized in all manuals. There may therefore be chapters within this document without any content. Chapters without content are identified as such.						
Styles used in the document	The following text styles may be used in the document: <table border="1"> <thead> <tr> <th>Style</th><th>Description</th></tr> </thead> <tbody> <tr> <td>Keys and buttons</td><td>Keys and buttons on the device are shown in bold type. Example: Example key.</td></tr> <tr> <td>➤ Instruction</td><td>➤ Instructions are indicated by an arrow ➤. Instructions must be followed. ➤ Example: ➤ Carry out instruction.</td></tr> </tbody> </table>	Style	Description	Keys and buttons	Keys and buttons on the device are shown in bold type . Example: Example key.	➤ Instruction	➤ Instructions are indicated by an arrow ➤. Instructions must be followed. ➤ Example: ➤ Carry out instruction.
Style	Description						
Keys and buttons	Keys and buttons on the device are shown in bold type . Example: Example key.						
➤ Instruction	➤ Instructions are indicated by an arrow ➤. Instructions must be followed. ➤ Example: ➤ Carry out instruction.						

Style	Description
1. Numbered instruction 2. ... 3. ...	Long passages containing instructions can be shown as numbered lists. Instructions must be followed. Example: 1. Carry out instruction.

Illustrations

The illustrations used in the documents may differ from the original if this does not have any influence on the function.

Importance of the instructions

These Instructions for Use are part of the accompanying documents and are an essential part of the device. They contain all the information necessary for operating the device.

The Instructions for Use must be carefully studied before attempting to operate the device.

Before the responsible organization may begin operating the device, the individual responsible for the operation must have been instructed by the manufacturer on how to use the device and must be thoroughly familiar with the contents of the Instructions for Use.

The device must only be operated by individuals who have been trained and certified in the proper operation and handling of the device.

Changes

Changes to technical documents will be released as new editions or supplements. In general, this manual is subject to change without notice.

Reproduction

Reproduction, even in part, is only permitted with written approval.

2.2 Significance of warnings



Warning

Advises the operator that failure to observe this information can result in personal injury.

2.3 Significance of notes



Note

Advises the operator that failure to observe this information can result in the following:

- Damage to the device
- Required functions will not be executed at all or not executed correctly

2.4 Significance of tips



Tip

Information providing useful tips for easy handling.

2.5 Brief description

Automated peritoneal dialysis (APD) is a treatment method normally performed overnight. Dialysis fluid is infused and drained by a so-called cycler.

The *sleep-safe harmony* technology

Important components of the *sleep-safe harmony* system are:

- Easy handling due to a clearly arranged color screen (touchscreen) for displaying and entering data.

- Automatic detection and connection of bags.
- Automatic inline flow heating of the dialysate to body temperature directly during the inflow phase.
- Options for profiling treatment parameters, for example, different glucose concentrations, inflow volumes and dwell durations can be programmed.
- Adapted APD: short dwell durations with low inflow volumes and long dwell durations with high inflow volumes
- Patient card: saving patient data, prescription data and treatment data on a mobile patient card.
Up to nine different prescriptions as well as the treatment reports of more than one year can be saved on the patient card.
- Patient connector with PIN technology reduces the risk of contamination during patient disconnection.

The *sleep•safe harmony* offers the possibility of performing CCPD, IPD, NIPD, Tidal Dialysis or PD-Plus therapies.

The device is classified as class IIb equipment (MDD).

2.6 Intended use

2.6.1 Intended purpose

The *sleep•safe harmony* has been designed for the treatment of patients suffering from renal insufficiency or terminal renal failure. It offers the possibility of performing CCPD, IPD, NIPD, Tidal Dialysis or PD-Plus therapies. The therapies available on the device include fluid removal and detoxification.

2.6.2 Specification of use

The device has been specified by the manufacturer for

- The treatment of patients with a body weight of 20 kg and more, irrespective of their age, under consideration of the specified technical data of the *sleep•safe harmony* (e.g., volume, delivery rates).
- Continuous treatments of up to 24 hours, with a typical overnight treatment duration of 7–10 hours.
- Operation in rooms that are suitable for peritoneal dialysis located in professional health care institutions, or for the home healthcare environment.

2.6.3 Side effects

Peritoneal dialysis therapy carries the risk of an inflammation of the peritoneum (peritonitis). In general, a lower incidence of peritonitis has been observed in APD therapy compared with CAPD therapy. Furthermore, infections of the catheter exit site may occur, which often progress to tunnel infections and, eventually, to peritonitis.

Distension and bloating (abdominal pain) may also occur in some patients. Shoulder pains and shortness of breath have also been observed, caused by an elevated diaphragm.

Depending on the dialysis solution used, a water-electrolyte imbalance can occur, e.g., potassium deficiency (hypokalemia). If too much fluid is removed during peritoneal dialysis, this may result in a decreased blood volume (hypovolemia) with attendant diminished blood pressure.

2.6.4 Contraindications

None known.

2.6.5 Interactions with other systems

None known.

2.6.6 Therapy restrictions

None known.

2.6.7 Target group



Warning

The device must only be installed, operated and used by individuals with the appropriate training, knowledge and experience and who are certified to have been trained.

2.7 Considerations for working on the device



Warning

Start-up, extensions, adjustments, calibrations, maintenance procedures, modifications or repairs must only be carried out by the manufacturer or persons authorized by the manufacturer.

All steps and information required for repairs to the device are included in the technical descriptions in the Service Manual.

For information on installation, refer to the appropriate chapter (see Chapter 9 on page 9-1).

For information on the Technical Safety Checks and maintenance procedures, refer to the appropriate chapter (see Chapter 11 on page 11-1).

Use only original spare parts.

For identifying and ordering spare parts, test equipment and tools, always use the electronic spare parts catalog.

Transport and storage (see Chapter 10 on page 10-1).

2.8 Expected service life

If the Technical Safety Checks are performed to the full extent specified and at the prescribed intervals, the safe operation of the device in the time between TSC procedures is guaranteed.

In addition, the manufacturer recommends that maintenance procedures be performed at the same time intervals to avoid device malfunctions caused by wear and tear.

With each Technical Safety Check, the "expected service life" according to IEC 60601-1 will therefore be prolonged until the next prescribed Technical Safety Check.

2.9 Duties of the responsible organization

The responsible organization has the following duties:

- Compliance with the national or local regulations concerning the installation, operation, use, and maintenance of the device.
- Compliance with the accident prevention regulations.
- Ensuring the proper and safe condition of the device.
- Ensuring the permanent availability of the Instructions for Use.

2.10 Operator responsibility



Warning

The following must be observed when entering parameters:

- The parameters entered must be verified by the operator, i.e., the operator must check that the values entered are correct.
- If the check reveals a deviation between the required parameters and the parameters displayed on the device, the setting must be corrected before activating the function.
- The actual values displayed must be compared with the prescribed target values.



Warning

The treatment parameters loaded from the *sleep•safe harmony* patient card must be checked by the operator for plausibility and compliance with the physician's prescription before being used.

For reasons of data protection, the operator is responsible for the safekeeping of the patient card.

The device must only be operated under the operating conditions specified by the manufacturer.

2.11 Disclaimer of liability



Warning

The device has been approved for use with certain consumables and accessories (see Chapter 8 on page 8-1). Should the responsible organization wish to use other consumables and accessories than those listed in the relevant chapter, their suitability must be verified beforehand by obtaining the appropriate information from the manufacturers.

The applicable legal regulations must be complied with.

The manufacturer does not assume any responsibility or liability for personal injury or other damage and excludes any warranty for damage to the device resulting from the use of non-approved or unsuitable consumables or accessories.

2.12 Warnings

2.12.1 Warnings, hygiene



Warning

Risk of infection

- Use aseptic technique when connecting the patient connector.
- You are recommended to wear a face mask, wash your hands and the spaces between your fingers with medical grade handwash and then apply a hand disinfection rub.
- Observe the hygiene practices of the dialysis center and the hygiene regulations in force.

The attending physician is responsible for proper patient training.

2.12.2 Warnings, therapy



Warning

Risk of injury through damaged line set.

Risk of infection

- When inserting the line set (*sleep•safe* Set and the connected solution bags), make sure there is no risk of damage by sharp-edged objects or by pets.



Warning

Patients suffering from severe ascites are critically ill patients who require medical supervision during treatment according to their particular clinical picture.

If these patients lose large amounts of ascitic fluid with their outflow, the appropriate therapeutic measures must be taken to counteract the possibility of a fall in blood pressure or hemodynamic instability.

2.12.3 Warnings, system



Warning

If the device is used outside the specified storage and operating conditions, the safe function of the device cannot be guaranteed.

This may result in damage to patients, operators or the environment.



Warning

Risk of injury through defective device.

Do not use the device, and terminate a therapy in progress, if any of the following defects are detected:

- Mechanical damage
- Defective power cable
- Other defects
- The device does not react as expected
- The performance characteristics deteriorate

➤ If the treatment is stopped due to an alarm (system error / device fault), follow the instructions of the attending physician.

➤ The device must be taken out of service.

➤ Inform the responsible organization or service support.



Warning

Risk of injury through swallowing small parts.

Choking hazard due to small parts.

➤ Keep loose small parts out of the reach of children.



Warning

Strangulation hazard.

➤ Make sure the cables and lines do not present a hazard for children.

2.12.4 Warnings, electrical



Warning

Risk of injury from electrical currents.

Electric shock

- The power cable must be laid so as to ensure that it cannot be damaged by sharp-edged objects or by pets.



Warning

- The use of power strips and extension cables is prohibited.



Warning

- Operators must not simultaneously touch the patient and the contacts of the plugs or sockets of the device. This can cause a patient leakage current.

2.12.5 Warnings, consumables and accessories



Warning

- The *sleep•safe* Set and the *sleep•safe* Set Plus are single use items. The sets are specially designed for the *sleep•safe harmony*. Only these sets must be used.



Warning

- Consumables must be discarded after treatment in compliance with the regulations for the disposal of potentially contaminated materials.

2.13 Addresses

Manufacturer Fresenius Medical Care AG & Co. KGaA
61346 Bad Homburg
Germany
Phone: +49 6172 609-0
www.fmc-ag.com


International service support Fresenius Medical Care
Deutschland GmbH
Service Support International
Hafenstrasse 9
97424 Schweinfurt
Germany
Phone: +49 9721 678-333 (Hotline)
Fax: +49 9721 678-130

Local service support

3 Design

3.1 Views

3.1.1 Front view



1 External key

This key is used to confirm the connection or disconnection of the patient as well as certain screen messages.

2 Status indicator

The indicator lights up red to indicate an alarm, and during the functional test.

The indicator lights up green to indicate correct operation.

3 Screen with integrated touch function

The screen displays treatment information and buttons that can be pressed.

4 On/Off key

This key is used to switch the *sleep•safe harmony* on or off.

5 Card slot

Card slot for the patient card.

The patient card contains the patient's individual prescriptions and treatment data.

6 Loading tray

The loading tray holds the *sleep•safe Set*.

7 Connector rail

The connectors of the solution bags are inserted into the connector rail.

3.1.2 Rear view

**1 Alarm output**

Can be used to connect an external alarm indicator.

2 LAN interface

Interface for data exchange. Here, optional device extensions can be connected.

3 Identification label

The identification label shows the *sleep•safe harmony* connection data.


4 Power cable

The power cable is used to connect the *sleep•safe harmony* to the power supply.

5 Power switch

The power switch is used to switch the electrical power on or off.

3.2 User interface

**A Status bar**

1 Help

2 Title of the action to be performed or status of the *sleep•safe harmony*

B Options panel

3 Tabs

4 Accept / Confirm button

C **Menu bar**

5 Device option menu

Depending on the treatment phase, you may not be able to select all of the device options.

6 Therapy option menu

7 Patient option menu

3.2.1 Screen colors


- **Color of control elements**

All the elements which can be selected (e.g., buttons, tabs) have a uniform operating concept.

blue: can be selected

green: active

gray: not active / cannot be selected



1 **Tabs**

green:

Tab selected.

The data for the selected tab will be displayed.


blue:

Tab not selected.

This tab can be selected.

3.3 General procedure for entering parameters

3.3.1 Entering text



➤ Use the keyboard displayed to enter the required text.

The other buttons have the following functions:


1. Shifts between upper case and lower case letters
2. Shifts between letters and digits
3. Space bar
4. Shows special characters
5. Deletes the last character

➤ When all the required parameters have been entered, press the button to confirm the data entered


or

➤ Press the button to discard the change.

3.3.2 Entering numbers




- Select the digit to be edited. The selected digit is shown in green.
- Use the and buttons to increase or decrease the digit.
- If necessary, select and edit the next digit.
- When all the required parameters have been entered, press the button to confirm the data entered or
- Press the button to discard the change.



A screen message with a green frame confirms that the data has been saved.

If the value entered is outside the valid range, the *sleep•safe harmony* will adjust this value to the next possible, permitted value.

- When all the required parameters have been entered, press the button to confirm the data entered or
- Press the button to discard the change.



If the value entered is outside the range which can be set by the *sleep•safe harmony*, the value entered will be rejected.


An appropriate screen message with a yellow frame will be displayed at the bottom of the screen.

The value entered must be within the valid range.

In this example, either the number of base cycles or the dwell duration can be reduced.

- Press the button to confirm the input
- or
- Press the button to discard the change.


3.4 Selecting and editing options or parameters



- Press the **Therapy** button.




- Use the and buttons to move the required option to the highlighted line in the center. Press on this line to select the option.
- Pressing the button will display the higher-level screen.




- Use the and buttons to move the required parameter to the highlighted line in the center. Press on this line to select the parameter.

The number of the selected list item (1) and the total number of list items (2) in this menu are displayed at the bottom of the screen on the right.

- Pressing the button will display the higher-level screen.



- Enter the required value (see Chapter 3.3.2 on page 3-7).
- Press the button to confirm the input
- or
- Press the button to discard the change.



A screen message with a green frame confirms that the data has been saved.

If the value entered is outside the valid range, the *sleep•safe harmony* will adjust this value to the next possible, permitted value.

- When all the required parameters have been entered, press the button to confirm the data entered
- or
- Press the button to discard the change.



The prescription data can be saved in two ways.

The prescription data is saved permanently (1) and can be used for future treatments

or

the prescription data is saved temporarily (2) and can only be used for the current treatment.

- Press the appropriate button.
- Press the  button to confirm the input
- or
- Press the  button to discard the change.

3.4.1 Selecting and editing the solution

There are two ways to select the solution:


- Selecting the solution directly from the list.
- Restricting the selection by selecting only a specific solution type, a specific glucose concentration, a specific calcium concentration (Ca^{++}) or a combination of these three selection options.



➤ Press the **Therapy** button.




➤ Use the and buttons to move the option **Edit prescription** to the highlighted line in the center. Press on this line to select the option.
➤ Pressing the button will display the higher-level screen.




- Use the and buttons to move the parameter **Solution** to the highlighted line in the center. Press on this line to select the parameter.
- Pressing the button will display the higher-level screen.

- **Selecting the solution directly**



- Use the and button to move the required solution to the highlighted line in the center. Press on this line to select the parameter.
- Press the button to discard the change.

- **Restricting the selection of solutions**




The buttons 1, 2 or 3 can be used to restrict the selection of the solutions.

1. Solution type
2. Glucose concentration
3. Calcium concentration


The procedure for restricting the selection is described using the glucose concentration as an example.

- Press button 2.
- Press the button to discard the change.



- Use the and buttons to move the required parameter to the highlighted line in the center. Press on this line to select the parameter.

- Pressing the button will display the higher-level screen.



The selected parameter is displayed in the appropriate field (1).

- In the restricted selection list, use the and buttons to move the required parameter to the highlighted line in the center. Press on this line to select the parameter.
- Press the button to discard the change.




Note

When restricting the selection of solutions in this way, certain parameter combinations may result in no matching solution being found.

In this case an appropriate message will be displayed.

3.5 Profiling

Some types of prescriptions permit parameter profiling. Here, different parameter values can be set for the individual treatment cycles.




➤ Press the **Therapy** button.



➤ Use the and buttons to move the option **Edit prescription** to the highlighted line in the center. Press on this line to select the parameter.

➤ Pressing the button will display the higher-level screen.



The icon (1) next to the parameters indicates which parameter can be profiled.

Parameter already profiled.


Parameter not yet profiled.

The parameters (2) are valid for a non-profiled prescription.

If a parameter is profiled, the respective parameter is not active and will be shown in gray.

➤ Select the **Details** tab.

➤ Pressing the button will display the higher-level screen.



The following three parameters can be profiled:


- Inflow volume
- Solution
- Dwell duration

➤ Press the appropriate button to profile the required parameter.

➤ Pressing the button will display the higher-level screen.




- Use the and buttons to move the required profile type to the highlighted line in the center. Press on this line to select the parameter.
- The button cancels the profiling.



Parameters that can be profiled are identified by the icon (1).

The treatment cycle for which the parameter value can be changed is highlighted (2).


- Use the and buttons to move the treatment cycle for which the parameter is to be edited to the center. Press on the highlighted parameter field to open the screen for entering or selecting the required value.



- Enter the required value (see Chapter 3.3.2 on page 3-7).
- Press the button to confirm the input

or

- Press the button to discard the change.



The changed value for this treatment cycle will be displayed.

A screen message with a green frame confirms that the data has been saved.


If the value entered is outside the valid range, the *sleep•safe harmony* will adjust this value to the next possible, permitted value.

An appropriate message will be displayed at the bottom of the screen.

- When all the required parameters have been entered, press the button to save the changes

or

- Press the button to discard the change.



The prescription data can be saved in two ways.

The prescription data is saved permanently (1) and can be used for future treatments

or

the prescription data is saved temporarily (2) and can only be used for the current treatment.

➤ Press the appropriate button.

➤ Press the button to confirm the input

or

➤ Press the button to discard the change.

4 Operation

4.1 Preparing the device

This chapter describes how to prepare the *sleep•safe harmony* for treatment.



Note

- Do not place solution bags on the *sleep•safe harmony*.



Note

- During outflow, the level difference between the *sleep•safe harmony* and the patient should not be significantly changed. A change between a lying and a sitting position does not have any negative effect.

4.1.1 Functional test



- Connect the *sleep•safe harmony* to the power supply.

- Switch the *sleep•safe harmony* on with the power switch.

After the power has been switched on, the device will first perform a functional test.


The  key will light up after approx. 1 minute.

- Press the  key as soon as it lights up.



Warning

- Only the patient card approved by Fresenius Medical Care and intended for use with the *sleep•safe harmony* may be used in conjunction with the *sleep•safe harmony*.
If the *sleep•safe harmony* is operated without a patient card, the treatment data stored in the *sleep•safe harmony* will be used for the treatment. On completion of the treatment, the treatment results will be saved in the *sleep•safe harmony*.
If the patient card was inserted into the *sleep•safe harmony* before the device was switched on, the treatment data stored on the patient card will be used for the next treatment. On completion of the treatment, the treatment results will be saved to the patient card.
- Do not remove the patient card from the *sleep•safe harmony* while treatment is in progress.



The *sleep-safe harmony* is ready for operation as soon as the screen shown on the left is displayed.

The following information will be displayed:

- The current date
- The software version
- The current time

➤ Press the nose of the smiley to move to the next operating step.

➤ If the *sleep-safe* Set is still in the *sleep-safe harmony* after a treatment has been prematurely terminated, press the button (1) on the right at the bottom of the screen to open the loading tray.

➤ For Service login, press the button (2) on the left at the bottom of the screen (only available for service support).

4.2 Preparing for treatment

4.2.1 Preparing the supplies



Warning

Risk of infection

- Use aseptic technique when connecting the patient connector.
- You are recommended to wear a face mask, wash your hands and the spaces between your fingers with medical grade handwash and then apply a hand disinfection rub.
- Observe the hygiene practices of the dialysis center and the hygiene regulations in force.
- Do not remove the single use items or consumables from their overwrap before this operating step is displayed by the *sleep•safe harmony*.
- Always make sure that the *sleep•safe* Set foil is not dented or damaged (deformed, torn, etc.).

The attending physician is responsible for proper patient training.

- Close windows and doors.
- Take off wrist watch and jewellery.
- Prepare the necessary supplies:
 1. *sleep•safe harmony*
 2. Patient card
 3. Solution bags in overwrap
 4. *sleep•safe* Set
 5. Drain line system
 6. Organizer with clip
 7. Disinfection cap

8. Face mask
9. Liquid soap (medical grade handwash)
10. Hand disinfectant



Note

➤ The solution bags used must have a temperature of 15–35 °C.



Warning

➤ Do not use the single use items or consumables if

- The overwrap is damaged.
- They are past their "use by" date
- The protective caps and closing caps have come off



➤ Check the following items on each solution bag:

- The name of the dialysis solution matches that in the prescription
- The glucose concentration and the calcium concentration match those specified in the prescription
- It is not yet past its "use by" date
- The overwrap is not damaged
- There is no leakage potential
- The dialysis solution is clear

- **Handling 5 liter double chamber bags**



- Before using the bags, check that both PEEL seams are intact.
- Open the overwrap and leave the bag on the lower part of the overwrap.
- Fold out the separating seam in the middle and the bag connector.
- Roll up the solution bag starting from the upper corner, which is diagonally opposite the bag connector, until the PEEL seam opens.



- Continue rolling up the solution bag until the PEEL seam of the small chamber opens completely.

All PEEL seams must now be open.

- Mix the dialysis solution and check that the bag is tight.

- **Handling 3 liter double chamber bags**



- Before using the bags, check that both PEEL seams are intact.
- Open the overwrap and leave the bag on the lower part of the overwrap.
- Roll up the solution bag starting at the upper corner until the PEEL seam opens.




- Then roll up the solution bag starting from the upper edge until the PEEL seam of the lower triangular section opens completely.

All PEEL seams must now be open.

- Mix the dialysis solution and check that the bag is tight.

4.2.2 Confirming the prescription



The **Prescription** tab shows the following information:

- Patient name
- Prescription name
- Maximum inflow volume
- Total treatment volume
- Expected treatment duration

➤ Select the **Solution** tab

or

➤ Press the  button to move to the next operating step and to start preparation.




The **Solution** tab shows the following information:

- Solution type
- Amount of glucose in %
- Amount of calcium in mmol/l
- Total solution volume in ml of the respective solution type

➤ Select the **Graphics** tab

or

➤ select the **Prescription** tab.



The **Graphics** tab shows a schematic of the prescription.

The vertical axis shows the volume in ml and the horizontal axis shows the time in minutes.

1. Permitted patient volume in ml (orange) (see Chapter 7.3.1.1 on page 7-9)
2. Permitted residual volume in ml (green) (see Chapter 7.3.1.2 on page 7-9)
3. Prescribed volume curve (blue)


➤ select the **Prescription** tab.

4.2.3 Inserting the *sleep•safe* Set




Warning

The treatment parameters loaded from the *sleep•safe harmony* patient card must be checked by the operator for plausibility and compliance with the physician's prescription before being used.



➤ Press the button to confirm the prescription.



The *sleep•safe harmony* will perform internal tests.

The loading tray opens.
During the test an audible alarm can be heard and the status indicator is red.
After about 5 seconds, first the connector rail, then the loading tray opens.

The next screen message will be displayed automatically.



Warning

- Do not use the *sleep•safe harmony* if the audible alarm is not heard or the status indicator is not lit during the internal test.
- Call service support.



Warning

- Do not place any objects on the *sleep•safe Set*.



- Open the overwrap of the *sleep-safe* Set and remove the *sleep-safe* Set.
- Place the *sleep-safe* Set in the open loading tray.
- Push the lines on the left and right into the line locators provided for this purpose.
- Place the lines in the line guides.
- The  button can be used to replay the screen animation.



Warning

- If the *sleep-safe* Set gets jammed when the loading tray closes, preventing the tray from being closed properly, the *sleep-safe* Set can be damaged.
- Use a new *sleep-safe* Set for the treatment.



- Press the  button to confirm that the *sleep-safe* Set is positioned correctly and to move to the next operating step.

4.2.4 Connecting the solution bags




Warning

➤ For reasons of hygiene, the solution bag connectors and the connector ends of the *sleep•safe* Set must only be touched on the outside.



Note

To prevent errors when checking the solution bag connectors, protect the loading tray slot from direct light.



The solution bags required for the treatment are shown on the screen. The color coding facilitates the correct selection of the solution bags.

- Take off the overwrap of the solution bags.
- If using a double chamber bag, mix the solution.
- Check the solution bags for leakage before use.
- Unroll the connecting line with the solution bag connector.
- Unscrew the closing cap of the solution bag connector.
- Place the connectors of the solution bags in a free tray port in front of the *sleep•safe* Set and push down.

Make sure all of the connectors are firmly secured in the tray ports.

- When using a 7.5 % polyglucose solution bag, handle the bag as shown in the operator training (see Chapter 4.2.4.2 on page 4-16).
- Pressing the  button will display the higher-level screen.
- The  button can be used to replay the screen animation.
- Press the  button to move to the next operating step.

Closing loading tray ...



The loading tray closes.

The next screen message will be displayed automatically.


Checking solution bag connectors ...



The volume and solution type of the inserted solution bags are compared with the prescription.

The next screen message will be displayed automatically.

4.2.4.1 Connected solution bags do not match the prescription



The loading tray will be opened if the connected solution bags do not match the prescription. The incorrect solution bags are identified by a red arrow.

- Check the solution bag and replace it as necessary.
- Press the  button to return to the previous operating step or
- press the  button to repeat the barcode reading.




The volume and solution type of the inserted solution bags are compared with the prescription.

The next screen message will be displayed automatically.

4.2.4.2 Continuing the connection of the solution bags

- **Confirming that double chamber bags have been mixed**



If a single chamber dialysis solution bag is used, this screen will not be displayed.


➤ Press the button to move to the next operating step.

- **Connecting a 7.5 % polyglucose solution bag**




Warning

➤ The *sleep•safe* Luer lock adapter must only be used in combination with a 7.5 % polyglucose solution bag, for the last inflow.



- When connecting a 7.5 % polyglucose solution bag, make sure you have the Luer lock adapter and the 7.5 % polyglucose solution bag ready to hand.
- Before breaking the cone of the 7.5 % polyglucose solution bag, make sure the Luer lock adapter is correctly connected.
- Press the button to move to the next operating step.



The connector rail will close when all the required connectors have been correctly inserted.

The solution bags will then be connected to the *sleep•safe* Set automatically.

4.2.5 Connecting the drain line system




Note

- The use of the *sleep•safe* drain line extension is recommended if 5 or 6 solution bags with 5 liters or 6 liters each are connected. The capacity of the drain line set is not sufficient for a treatment with 6 solution bags.


**Note**

- The outlet of the drain line must never be more than 2 meters below or above the *sleep•safe harmony*.



- After connecting the solution bags, connect the drain line.
- Unroll the drain line on the right of the *sleep•safe* Set.
- Connect the yellow connector of the *sleep•safe* Set to the drain line or the drain line set.
- When using the drain line set, make sure that the clamp of the drain line is open and that the clamps of the emptying lines on the drain line set are closed.
- The button can be used to replay the screen animation.
- Press the button to move to the next operating step.


4.2.6 Priming



The **sleep-safe** Set, the patient line and the drain line will be primed automatically.

The next screen message will be displayed automatically.

4.2.6.1 Checking the fill level in the patient line




- If priming has been completed, but the patient line is not yet fully primed, it can be filled up with dialysis solution by pressing the  button.
- Press the  button to move to the next operating step.

4.2.7 Confirming treatment data



Warning

➤ The operator must check the treatment data (maximum inflow volume, treatment volume and treatment duration) for plausibility before starting the treatment.



Check the treatment data for plausibility.



➤ Press the key to confirm the treatment data.

4.3 Starting the treatment

4.3.1 Connecting the patient



Warning

Risk of infection


- Use aseptic technique when connecting the patient connector.
- You are recommended to wear a face mask, wash your hands and the spaces between your fingers with medical grade handwash and then apply a hand disinfection rub.
- Observe the hygiene practices of the dialysis center and the hygiene regulations in force.

The attending physician is responsible for proper patient training.



Warning

- Treatment should be performed immediately after priming is complete.



- Once priming is complete, the patient should be connected immediately.

The patient's name will be displayed.

- Connect the patient line as described below.
- The button can be used to replay the screen animation.



- Unroll the patient line on the left of the *sleep-safe* Set.
- Insert the patient connector in the clip.
- Place the catheter extension system connector into the right holder of the organizer (if left-handed, place in the left holder).



- Put on the face mask.
- Disinfect your hands and dry them carefully.



- Unscrew and discard the protective cap of the *sleep-safe* patient connector.



- Unscrew the catheter extension system connector from the disinfection cap.




- Screw the catheter extension system connector directly onto the patient connector of the *sleep•safe* Set.
- Open the white clamp on the catheter extension.



Warning

- Only the patient whose name is displayed on the screen must be connected to the *sleep•safe harmony*.



➤ Press the button to confirm the connection and to start the treatment.

4.3.2 Starting the treatment



Treatment is started.

The treatment will start with an initial outflow.

4.3.3 Terminating the initial outflow




Warning

After the initial outflow has been completed normally and the alarm has been confirmed by the operator, the *sleep•safe harmony* assumes that the draining of the peritoneal cavity has been completed. Volume monitoring will be restarted.




If an initial outflow has been set for the treatment, this message will be displayed when the device detects that the draining of the peritoneal cavity has been completed.

- Press the  button to silence the audible signal.



The initial outflow volume achieved and the initial outflow target volume will be displayed.

- Press the  button to repeat the initial outflow.
- Press the  button only if you wish to end the initial outflow and start the subsequent inflow phase.




This warning appears when you end the initial outflow.

- Do not press the key unless there is definitely no fluid left in the peritoneal cavity.
The first inflow will be started.
- Press the button to continue the initial outflow.



4.3.4 Treatment overview



During the treatment, the following parameters will be displayed on the treatment overview screen:

1. Current time
2. Treatment progress in %
3. Expected end of treatment
4. Current total volume balance
5. Volume in patient
6. Progress of the current treatment phase
7. Treatment cycle

Negative values describe a weight loss (ultrafiltrate has been generated). Positive values describe a weight gain (resorption).

4.4 Terminating the treatment

4.4.1 Disconnecting the patient




Warning

Risk of infection

- Use aseptic technique when disconnecting the patient connector.
- Wearing a face mask and hand disinfection is recommended.
- Observe the hygiene practices of the dialysis center and the hygiene regulations in force.


The attending physician is responsible for proper patient training.



Perform the steps for disconnection as detailed below.

1. Release the PIN
2. Close the clamp
3. Disconnection
4. Confirmation

➤ The  button can be used to replay the screen animation.



- Turn the blue knob on the patient connector one quarter of a turn clockwise.



- Push the blue knob on the patient connector fully in.

This will automatically close the catheter extension with the PIN.

- Close the white clamp on the catheter extension.



- Open the overwrap of the new disinfection cap.
- Place the new disinfection cap in the left holder of the organizer (if left-handed, place the disinfection cap in the right holder).
- Insert the patient connector into the clip of the organizer.
- Unscrew and discard the closing cap of the new disinfection cap.




- Unscrew the catheter extension from the patient connector of the *sleep•safe* Set and screw the catheter extension with the PIN immediately onto the new disinfection cap.



- Pull the closed catheter extension straight (without turning it) out of the organizer.




- Press the  key to confirm the disconnection.



The treatment results will be saved to the patient card or in the device. If the patient card needs to be removed, e.g., for visiting a physician, this must only be done after the device has been switched off.

4.4.2 Treatment results



Select the **Therapy** menu to display the treatment report or press the  button to terminate the treatment.

The results of the last treatment will be displayed.

1. Total of all inflow volumes plus last inflow
2. Volume balance without last inflow and without initial outflow
Negative values describe a weight loss (ultrafiltrate has been generated). Positive values describe a weight gain (resorption).
3. Total treatment duration
4. Number of interruptions
5. Reduction of the dwell duration of the total treatment (as reduced automatically within the permitted limits) (see Chapter 7.3.2 on page 7-10).
6. Reduction of the inflow volume of the total treatment (as reduced automatically within the permitted limits) (see Chapter 7.3.1.3 on page 7-10).

➤ Press the  button to exit the **Treatment complete** screen and to move to **Draining**.




Warning

A repeated reduction of the dwell duration or of the inflow volume may result in an inadequate treatment.

In this case, please inform the physician.




- Select the **Treatment report** option.
- Pressing the  button will display the higher-level screen.



The **Details** tab shows the results of the treatment cycles and the total treatment results.

1. Treatment cycle
2. Inflow volume in ml and inflow duration in minutes
3. Patient volume in ml and dwell duration in minutes
4. Outflow volume in ml and outflow duration in minutes
5. Solutions
6. Total treatment duration in hours and minutes
7. Balance in ml
8. Total treatment volume in liters

➤ Select the **Interruptions** tab
or
➤ Pressing the  button will display the higher-level screen.




The **Interruptions** tab shows all the cautions and alarms, the time that they occurred, and the treatment phase during which they occurred.

1. Treatment phase and cycle
2. Number and description of the interruption
3. Operator response to the interruption
4. Time when the interruption occurred

The significance of the interruptions is explained in a separate chapter (see Chapter 5 on page 5-1).

- Select the **Graphics** tab
- or
- Pressing the button will display the higher-level screen.



The **Graphics** tab shows a schematic of the treatment results.

The vertical axis shows the volume in ml and the horizontal axis shows the time in minutes.

1. Permitted patient volume in ml (orange) (see Chapter 7.3.1.1 on page 7-9)
2. Permitted residual volume in ml (green) (see Chapter 7.3.1.2 on page 7-9)
3. Prescribed volume curve (blue)
4. Actual volume curve (red)
5. Volume removed curve (yellow) (Example as shown: -1200 ml indicates a total removed volume of 1200 ml).


➤ Select the **Interruptions** or the **Details** tab.

4.4.3 Draining the line set



Note

➤ Do not place the solution bags above the level of the *sleep•safe harmony*, if draining has been prematurely stopped or has been skipped and there is still fluid left in the solution bags.



The solution bags and the line set are drained automatically.

The *sleep-safe harmony* prepares for unloading the *sleep-safe Set*.

➤ The  button can be used to skip draining.

The loading tray will open automatically when the automatic draining of the line set has been completed, or if draining was stopped manually.

While the loading tray opens, a tone sequence can be heard.

4.4.4 Removing the *sleep-safe Set*



Note

➤ Do not pull on the connector line of the solution bags as this may damage the *sleep-safe Set* and cause a contamination of the *sleep-safe harmony* with dialysis solution.



➤ Remove the *sleep-safe Set* from the loading tray. Pull the patient line and the drain line simultaneously upwards until the lines come off the line holders.

➤ The  button can be used to replay the screen animation.

➤ Press the  button to confirm the removal of the set.



Warning

- Consumables must be discarded after treatment in compliance with the regulations for the disposal of potentially contaminated materials.

Closing loading tray ...



The loading tray closes.

The *sleep•safe harmony* is checked and switches to standby mode.

If necessary, the *sleep•safe harmony* can be switched off with the power switch.

- Dispose of the *sleep•safe* Set, the solution bags and the drain line system.

Device can be turned off



Goodbye




Therapy




4.5 Special operational features of the *sleep-safe harmony*

4.5.1 Personalizing the *sleep-safe harmony* for a patient


When the *sleep-safe harmony* is used for the first time for a patient, the device or the patient card must be configured for the patient.




➤ Press the button.




➤ Enter the password.
 ➤ Press the button to confirm the input and to move to the next operating step.
 ➤ Pressing the button will display the higher-level screen.




- Press the **Select language** button.
- Select the required language.
- Press the button to confirm the input and to move to the next operating step.
- Press the button to return to the previous operating step.




- Press the **Enter first name** button.
- Enter the first name.
- Press the button to confirm the input and to move to the next operating step.
- Press the button to return to the previous operating step.




- Press the **Enter last name** button.
- Enter the last name.
- Press the  button to confirm the input and to move to the next operating step.
- Press the  button to return to the previous operating step.




- Press the **Enter date of birth** button.
- Enter the date of birth.
- Press the  button to confirm the input and to move to the next operating step.
- Press the  button to return to the previous operating step.




- Press the **Select prescription** button.
- Press the  button to confirm the input and to move to the next operating step.
- Press the  button to return to the previous operating step.




- Select the required prescription type and enter a name for this prescription.
- Pressing the  button will display the higher-level screen.



- Press the button to confirm the input and to move to the next operating step.
- Press the button to return to the previous operating step.



- Check all data entered.
- Once the data has been confirmed, the personalization is complete and the patient card is assigned to the patient. The first name, surname and date of birth can then no longer be edited on the *sleep•safe harmony*.
- Press the button to create a patient card.
- Press the button to return to the previous operating step.



➤ Check and edit the patient parameters and the prescription data, if necessary.

Editing the patient parameters (see Chapter 4.6.1.3 on page 4-54).


Editing the prescription data (see Chapter 4.6.2.3 on page 4-65).

➤ Press the button to create a patient card.



➤ On completion of the treatment or after switching the *sleep•safe harmony* off, the patient data and, depending on the selection, also the prescription data, will be saved to the patient card.

4.5.2 Preparing the *sleep-safe harmony* for a treatment without a patient card



➤ Press the button to start the treatment with a patient card

or

➤ Press the button to start the treatment without a patient card.

➤ The button can be used to replay the screen animation.



Warning

If the treatment was started without a patient card, the patient name must be confirmed.

If the patient name displayed is not the name of the patient to be treated, start personalization.




➤ Press the button only if you need to personalize the *sleep-safe harmony*.

This will create a patient card with standard values (see Chapter 4.5.1 on page 4-39).

➤ Press the button to confirm the patient name displayed.

4.5.3 Removing the *sleep-safe* Set after premature termination of a treatment



➤ Press the button (1) on the right at the bottom of the screen to open the loading tray.



➤ Press the button to confirm the opening of the loading tray.



Note

➤ Do not pull on the connector line of the solution bags as this may damage the *sleep-safe* Set and cause a contamination of the *sleep-safe harmony* with dialysis solution.



- Remove the *sleep•safe* Set from the loading tray. Pull the patient line and the drain line simultaneously upwards until the lines come off the line holders.
- Press the button to close the loading tray and to switch the *sleep•safe harmony* to standby mode.



Warning

- Consumables must be discarded after treatment in compliance with the regulations for the disposal of potentially contaminated materials.

4.5.4 Disconnection after premature termination of a treatment




Warning

Risk of infection

- Use aseptic technique when disconnecting the patient connector.
- Wearing a face mask and hand disinfection is recommended.
- Observe the hygiene practices of the dialysis center and the hygiene regulations in force.

The attending physician is responsible for proper patient training.



This screen message will only be displayed if treatment has been prematurely terminated and the patient has not been disconnected.

Perform the steps for disconnection

1. Release the PIN
2. Close the clamp
3. Disconnection
4. Confirmation

as detailed below.

➤ The  button can be used to replay the screen animation.

➤ Turn the blue knob on the patient connector one quarter of a turn clockwise.





- Push the blue knob on the patient connector fully in.

The catheter extension will automatically be closed by the PIN.

- Close the white clamp on the catheter extension.



- Open the overwrap of the new disinfection cap.
- Place the new disinfection cap in the left holder of the organizer (if left-handed, place the disinfection cap in the right holder).
- Insert the patient connector into the clip of the organizer.
- Unscrew and discard the closing cap of the new disinfection cap.



- Unscrew the catheter extension from the patient connector of the *sleep•safe* Set and screw the catheter extension with the PIN immediately onto the new disinfection cap.



➤ Pull the closed catheter extension straight (without turning it) out of the organizer.



➤ Press the  key to confirm the disconnection.

4.5.4.1 Switching the *sleep-safe harmony* off




➤ Press the  key.

The *sleep-safe harmony* will switch off after 5 seconds.

4.6 Options / changing data prior to treatment

4.6.1 Patient options

4.6.1.1 Selecting the language



➤ Select the  menu.



➤ Select **Language**.


➤ Pressing the  button will display the higher-level screen.




- Select the required language.
- Pressing the button will display the higher-level screen.

4.6.1.2 Editing personal data

This function can only be accessed by the clinical staff.



- Select the menu.



➤ Select **Personal data**.

➤ Pressing the  button will display the higher-level screen.



The following **Personal data** will be displayed:

- First name
- Last name
- Date of birth

The following **Personal data** can be edited:

- Gender
- Weight

➤ Pressing the  button will display the higher-level screen.

4.6.1.3 Editing patient parameters

This function can only be accessed by the clinical staff.

If patient parameters are restricted, the prescriptions will automatically be adapted to the new patient parameters.



Warning

➤ After editing, the operator must check the treatment data (permitted patient volume, treatment volume and treatment duration) for plausibility.

The preselected parameters (e.g., preselected inflow volume) will be used when creating a new prescription.



Warning

The maximum parameters (e.g., maximum inflow volume) represent the limits for all prescriptions.


A change of the maximum parameters may influence existing prescriptions.

When changing the maximum parameters, these must be checked for compliance with the prescription.

Patient parameters can only be edited by the attending physician.

The following patient parameters can be set:

- Preselected inflow volume
- Maximum inflow volume
- Preselected inflow rate
- Maximum inflow rate
- Preselected outflow rate
- Maximum outflow rate
- Preselected dwell duration
- Maximum dwell duration
- Permitted patient volume
(see Chapter 7.3.1.1 on page 7-9)
- Permitted residual volume
(see Chapter 7.3.1.2 on page 7-9)
- Permitted dwell duration reduction
(see Chapter 7.3.2.1 on page 7-10)
- Permitted inflow volume reduction
(see Chapter 7.3.1.3 on page 7-10)
- Catheter performance
(see Chapter 7.3.2.2 on page 7-11)
- Set type (function cannot be selected)
- Additional outflow
- Access level




➤ Select the  menu.



➤ Select **Patient parameters**.


➤ Pressing the  button will display the higher-level screen.



The preselected and maximum **Patient parameters** can be selected and edited.

➤ Pressing the  button will display the higher-level screen.

4.6.1.4 Changing the additional outflow




➤ Select the  menu.




➤ Select **Patient parameters**.

➤ Pressing the  button will display the higher-level screen.



➤ Select **Additional outflow**.

➤ Pressing the  button will display the higher-level screen.




➤ The additional outflow option can be set as follows:

- No
- Yes (with sound)
- Yes (without sound)


Selecting the option Yes (without sound) initiates a treatment pause. Confirming the treatment pause continues the treatment.

➤ Pressing the  button will display the higher-level screen.

4.6.1.5 Access levels



➤ Select the  menu.




➤ Select **Patient parameters**.

➤ Pressing the button will display the higher-level screen.

➤ Select **Access level**.


➤ Pressing the button will display the higher-level screen.



- The following **access levels** can be selected:
 - No changes
 - Select prescription
 - Edit prescription
- Pressing the  button will display the higher-level screen.

4.6.2 Therapy options

4.6.2.1 Creating a prescription




- Select the **Therapy** menu.



➤ Select **Create prescription**.


A maximum of nine prescriptions can be stored in the *sleep•safe harmony*. If this number is exceeded, an appropriate screen message will be displayed.

➤ Pressing the  button will display the higher-level screen.




➤ Select the prescription to be created.

➤ Pressing the  button will display the higher-level screen.



- Enter the name for the prescription.
- Press the  button to confirm the change.
- Press the  button to discard the change.



- Enter the parameters for the new prescription you have created.


4.6.2.2 Selecting a prescription



➤ Select the **Therapy** menu.




➤ Press **Select prescription** to display the prescriptions available for selection.
➤ Pressing the  button will display the higher-level screen.



- Select the required prescription.
- Pressing the  button will display the higher-level screen.

4.6.2.3 Editing a prescription




- Select the **Therapy** menu.




➤ Select **Edit prescription**.

➤ Pressing the  button will display the higher-level screen.



➤ Enter the parameters for the prescription to be edited.

4.6.2.4 Deleting a prescription




➤ Select the **Therapy** menu.



➤ Select **Delete prescription**.

➤ Pressing the  button will display the higher-level screen.




- Select the prescription to be deleted.
- Press the selected prescription to delete this prescription.
- Pressing the  button will display the higher-level screen.

4.6.2.5 Displaying treatment reports



- Select **Show treatment reports**.
- Pressing the  button will display the higher-level screen.




Treatment results					35
	1	2	3	4	5
	14.03.2015	14.43 l	-1108 ml	08h 42min	2
	13.03.2015	13.98 l	-1128 ml	09h 01min	1
	12.03.2015	14.13 l	-1185 ml	09h 06min	3
	11.03.2015	14.28 l	-1023 ml	08h 46min	4
	10.03.2015	13.59 l	-894 ml	09h 02min	3

A list of the last treatments is displayed.

1. Start of treatment date
2. Total of all inflow volumes plus last inflow
3. Volume balance
If treatment is terminated prematurely, this value does not show the actual volume balance of the patient, but only the volume removed within the treatment time. (Example as shown: -1000 ml indicates a total removed volume of 1000 ml).
4. Total treatment duration
5. Number of interruptions

➤ Select a treatment to view further details.

➤ Pressing the  button will display the higher-level screen.




Treatment report						138
	Details	Interruptions	Graphics			
	#	ml min	ml min	ml min	ml min	
	1	2001 ml 12 min	2001 ml 60 min	2058 ml 19 min	balance 1.50% Ca++ 1.75 mmol/l	
	2	1998 ml 13 min	2000 ml 60 min	2095 ml 19 min	balance 1.50% Ca++ 1.75 mmol/l	
	3	2001 ml 13 min	2001 ml 60 min	2182 ml 22 min	balance 1.50% Ca++ 1.75 mmol/l	
		1500 ml 9 min			balance 2.30% Ca++ 1.75 mmol/l	
		7.50 l	-345 ml	04 h 59 min		

➤ Details of a selected treatment report (see Chapter 4.4.2 on page 4-32).

➤ Pressing the  button will display the higher-level screen.

4.6.2.6 Editing the prescription schedule




➤ Select the **Therapy** menu.



➤ Select **Edit prescription schedule**.

➤ Pressing the  button will display the higher-level screen.




The treatment days of the prescription in the highlighted line in the center (3) are indicated by a green checkmark in the weekly overview (1).

If the standard prescription is selected, this is indicated in the Default field (2) by a green checkmark.

If no prescription is selected for a day of the week, the standard prescription will automatically be used.

- Select the prescription for which the weekdays are to be defined.
- Pressing the button will display the higher-level screen.




➤ Select the days of the week on which the selected prescription is to be performed.

- If the selected prescription is to be used as the standard prescription, also select the **Default** check box.
- Press the button to confirm the data entered.
- Press the button to discard the change.

4.6.3 Device options

4.6.3.1 Setting the date and time

Date and time can only be set when no treatment is being performed.




➤ Select the  menu.




➤ Select **Date and time**.

➤ Pressing the  button will display the higher-level screen.



- Select **Date, Time or Time format** and set the date, the time or the time format.
- Press the  button to confirm the change.
- Pressing the  button will display the higher-level screen.

4.6.3.2 LAN setup data




- Select the  menu.



➤ Select **LAN setup data**.

➤ Pressing the button will display the higher-level screen.



➤ Select **IP address, Netmask and Device name** and enter the data.


➤ Pressing the button will display the higher-level screen.

4.6.3.3 Relocating the *sleep•safe harmony* within a building

During preparation, after the line set has been completely primed, the device can be set to relocation mode, to move the device to a different location. This allows the preparation of the device and the treatment of the patient to take place in different rooms.

**Note**

The **Relocate device** function is only possible with a **sleep•safe Set Plus**.



➤ Select the



➤ Select **Relocate device**.

➤ Pressing the




- Press the button to confirm relocation.
- Press the button to return to the previous operating step.




Warning


- The patient line must be sealed using aseptic technique.



- Seal the patient line.
- Press the button to confirm that the patient line is sealed.
- Press the button to return to the previous operating step.




- Disconnect the *sleep-safe harmony* from the power supply.



- The *sleep-safe harmony* can now be moved to a different location.
- After a maximum of 15 minutes, the *sleep-safe harmony* must be reconnected to the power supply.

The remaining time is displayed on the bottom of the screen.



If the *sleep-safe harmony* has not been reconnected to the power supply, the screen on the left will be displayed.


➤ The *sleep-safe harmony* must now be connected to the power supply immediately, as the *sleep-safe harmony* will otherwise shut down.



Preparation for treatment can be continued.

➤ Press the button to confirm the patient name and the prescription.

4.6.3.4 Setting the volume




➤ Select the  menu.



➤ Select **Volume**.

➤ Pressing the  button will display the higher-level screen.



- Select the required volume for all the sounds (except audible alarms).
- Press the **Test** button to test the sound volume you have set.
- If the **Key sound activated** option is selected, this is indicated by a checkmark. The device will then emit a tone each time a key or button is pressed.
- Press the button to confirm the input


or

- Press the button to discard the change.

4.6.3.5 Setting the brightness



- Select the menu.



- Select **Brightness**.
- Pressing the  button will display the higher-level screen.



Set the required brightness for the display and all indicator lights (except for the red light of the **status indicator**).

- Press the  button to confirm the input
- or
- Press the  button to discard the change.

4.6.3.6 Setting the screen saver




➤ Select the  menu.




➤ Select **Screen saver**.

➤ Pressing the  button will display the higher-level screen.



- Set the screen saver timeout ("Waiting time").
- If the **Screen saver activated** option is selected, this is indicated by a checkmark. The screen saver function is then active.
- Press the button to confirm the input
- or
- Press the button to discard the change.

4.6.3.7 Displaying device information



- Select the menu.



➤ Select **Device information**.


➤ Pressing the  button will display the higher-level screen.

Device information			
		Pump 1	Pump 2
Pressure sensor	mbar	27	-59
Static patient pressure	mbar	-79	-77
Pressure sensor offset	mbar	-10	0
Maximum position	µm	31648	30897
Minimum position	µm	13337	13095
Specified flow rate	ml/min	230	
Air detection	%	40.1	
Inflow volume	ml	0.0 - 2000.0	
Outflow volume	ml	34.4 - 2147483.6	

Device information is displayed over several pages.


➤ Pressing the  button will display the next page.

➤ Pressing the  button will display the higher-level screen.



Device information is displayed over several pages.

- Pressing the button will display the next page.
- Pressing the button will display the previous page.
- Pressing the button will display the higher-level screen.



Device information is displayed over several pages.

- Pressing the button will display the previous page.
- Pressing the button will display the higher-level screen.

4.6.3.8 Clinic login

This function can only be accessed by the clinical staff.




➤ Select the  menu.



➤ Select **Clinic login**.


➤ Pressing the  button will display the higher-level screen.



- Enter the password.
- Press the  button to confirm the input and to move to the next operating step.
- Pressing the  button will display the higher-level screen.

4.6.3.9 Service login

This function can only be accessed by service support.



- Select the  menu.



➤ Select **Service login**.

➤ Pressing the  button will display the higher-level screen.




➤ Enter the password.


➤ Press the  button to confirm the input and to move to the next operating step.

➤ Pressing the  button will display the higher-level screen.

4.6.3.10 Clinic / Service logout




➤ Select the  menu.




➤ Select **Log out**.

➤ Pressing the  button will display the higher-level screen.



- Press the button to log out the clinical or service support staff currently logged in.
- Press the button to discard the change.

4.7 Therapy options during treatment



There are two ways to select the therapy options.


➤ Select the **Therapy** menu or press the  button.

The following therapy options are available.

- Pause treatment
- Pause with disconnection
- Manual outflow
- Skip phase
- Terminate treatment
- Treatment report




4.7.1 Pause treatment




➤ Select **Treatment options**.

➤ Pressing the  button will display the higher-level screen.




➤ Select **Pause treatment** to interrupt the treatment (pause) for a short period of time.

➤ Pressing the  button will display the higher-level screen.



- Press the button to confirm the treatment pause.
- Pressing the button will display the previous page.



The duration of the pause will be displayed.


- Press the button to continue the treatment.

4.7.2 Pause with disconnection



➤ Select **Treatment options**.


➤ Pressing the  button will display the higher-level screen.



➤ Select **Pause with disconnection** to interrupt the treatment and to disconnect the patient line.

➤ After the pause, reconnect the patient line.

➤ Pressing the  button will display the higher-level screen.




- Press the button to confirm the treatment pause.
- Pressing the button will display the previous page.

- **Disconnecting the patient line**



Warning

- The patient line must be sealed using aseptic technique.



- Disconnect the patient line as shown in the screen animation (see Chapter 4.4.1 on page 4-28).
- The button can be used to replay the screen animation.

- **Connecting the patient line**



The duration of the pause will be displayed.

➤ Press the button to continue the treatment.




Warning

Risk of infection


- Use aseptic technique when connecting the patient connector.
- You are recommended to wear a face mask, wash your hands and the spaces between your fingers with medical grade handwash and then apply a hand disinfection rub.
- Observe the hygiene practices of the dialysis center and the hygiene regulations in force.

The attending physician is responsible for proper patient training.




- Connect the patient line as shown in the screen animation (see Chapter 4.3.1 on page 4-21).
- The button can be used to replay the screen animation.

4.7.3 Manual outflow




- Select **Treatment options**.
- Pressing the button will display the higher-level screen.



- Select **Manual outflow** to initiate a manual outflow.
- Pressing the  button will display the higher-level screen.



- Press the  button to start the manual outflow.
- Pressing the  button will display the higher-level screen.




The manual outflow will be performed.

- The button continues the treatment.




The manual outflow is interrupted.

- Touch the screen to confirm the interruption.



When the manual outflow is complete, the following screen will be displayed.

- Press the button to continue the manual outflow.
- Press the button only if you wish to end the manual outflow.




The manual outflow is complete.

- Press the **Continue treatment** button or the **Terminate treatment** button (see Chapter 4.4.1 on page 4-28).


4.7.4 Skip phase



- Select **Treatment options**.
- Pressing the button will display the higher-level screen.



- Select **Skip phase** to skip the current treatment phase.
- Pressing the button will display the higher-level screen.



- Press the button to skip the current treatment phase and to start the next treatment phase.
- Pressing the button will display the higher-level screen.

4.7.5 Terminating the treatment




- Select **Treatment options**.
- Pressing the button will display the higher-level screen.



- Select **Terminate treatment** to terminate the current treatment immediately and to disconnect the patient line.
- Pressing the  button will display the higher-level screen.



- Press the  button to terminate the treatment.
- Pressing the  button will display the higher-level screen.



➤ Perform the steps for disconnection

1. Release the PIN
2. Close the clamp
3. Disconnection
4. Confirmation

as detailed below (see Chapter 4.4.1 on page 4-28).

➤ The  button can be used to replay the screen animation.

4.7.6 Displaying the treatment report



➤ Select **Treatment report** (see Chapter 4.4.2 on page 4-32).


➤ Pressing the  button will display the higher-level screen.

4.8 Performing a PD-Plus treatment

4.8.1 Connecting the patient

Connect the patient (see Chapter 4.3.1 on page 4-21).

4.8.2 Starting the PD-Plus cycle



During the treatment, the following parameters will be displayed on the treatment overview screen:

1. Current time
2. Treatment progress in %
3. Expected end of treatment
4. Current total volume balance
5. Volume in patient
6. Progress of the current treatment phase
7. Treatment cycle

4.8.2.1 Stopping the initial outflow




Warning

After the initial outflow has been completed normally and the alarm has been confirmed by the operator, the *sleep•safe harmony* assumes that the draining of the peritoneal cavity has been completed. Volume monitoring will be restarted.




If an initial outflow has been set for the treatment, this message will be displayed when the device detects that the draining of the peritoneal cavity has been completed.

- Press the  button to silence the audible signal.



The outflow volume achieved (1) and the target volume (2) will be displayed.


- Press the  button to repeat the initial outflow.
- Press the  button only if you wish to end the initial outflow.



- Press the  button to repeat the initial outflow.



- Do not press the  key unless there is definitely no fluid left in the peritoneal cavity.



The PD-Plus inflow is followed by the PD-Plus dwell duration during which the *sleep•safe harmony* treatment is paused.

- Disconnect the patient (see Chapter 4.4.1 on page 4-28).
- Press the  button to terminate the treatment pause.



Warning


In the time before treatment is started, make sure that the temperature of the fluid in the patient line does not rise above 40 °C due to external influences (exposure to sun, radiators, etc.).

4.8.3 Connecting the patient for nighttime treatment

Connect the patient (see Chapter 4.3.1 on page 4-21).

5 Alarm processing

5.1 Screen layout



- 1 **Help button** (function currently not available.)
- 2 **Screen message number**
The screen message number permits an unambiguous analysis of the error and provides the contact in the dialysis center or the local service support organization with additional troubleshooting information.
- 3 **Screen message title**
- 4 **Information text**
- 5 **Options panel**
The following buttons may be displayed on an alarm screen:

Confirm button 
 Button for confirming the screen message.

Repeat button 

Button to repeat the current test step.

Skip phase button 

Button for terminating the current treatment phase and switching to the next treatment phase.

Terminate treatment button 

Button for terminating the treatment.

5.2 Resetting the audible alarm



If an audible signal or an audible alarm is sounded, it can be silenced for 6 minutes by touching the screen.




➤ If an audible signal or an audible alarm is sounded and the "Reset audible alarm" screen is not displayed, it can be silenced by pressing the  key.

5.3 Color identification of the screen messages


The messages are organized in three categories which can be identified by their frame color.

**Yellow frame,
information (example)**



- Information and treatment-related assistance which the operator must always observe and follow as necessary.
- Treatment can be continued.

**Orange frame, caution
(example)**




- The title bar contains a brief description of the cause of the problem. Follow the steps displayed on the screen to correct the problem.
- Treatment can be continued.

**Red frame, alarm
(system error / device
fault)**

**Red frame, device
fault (example)**



- Faults which only affect the *sleep•safe harmony*.
- Depending on the alarm situation, treatment can be continued or will be stopped.

Red frame, system error (example)

- Errors which affect the entire system, including lines, solution bags and drainage.
- Depending on the alarm situation, treatment can be continued or will be stopped.

5.4 Screen messages

Message	Cause	Action required
A001	Internal error	<ul style="list-style-type: none"> – Disconnect the patient. – Terminate and restart the treatment, if necessary.
12T 55T 110	Heater error	<ul style="list-style-type: none"> – The <i>sleep-safe harmony</i> should not be exposed to direct sunlight. – The ambient temperature is too high. – Switch the <i>sleep-safe harmony</i> off with the power switch, allow the device to cool down for several minutes, then switch it on again.
61T	Error during the initial test	<ul style="list-style-type: none"> – Remove the <i>sleep-safe Set</i> by pressing the button on the right at the bottom of the screen (see Chapter 4.5.3 on page 4-46). Restart the <i>sleep-safe harmony</i>.
62T	Error during the initial test	<ul style="list-style-type: none"> – Device is running outside specified operating conditions, or time for adjustment to ambient temperature not observed (see Chapter 9.2 on page 9-2).
72T	Error during the initial test	<ul style="list-style-type: none"> – Check the patient line and the catheter extension for kinks and closed clamps. – Mix the solution in the double chamber bags. – Air detected in solution bag: Suspend the solution bag with the solution bag line pointing down and continue the treatment.
93	Pressure in patient line too high	<ul style="list-style-type: none"> – Check the patient line and the catheter extension for kinks and closed clamps. – Observe the patient position (see Chapter 9.3 on page 9-2). – Press the  button to terminate the current treatment phase and to switch to the next treatment phase.

Message	Cause	Action required
93T	Patient position monitoring	<ul style="list-style-type: none"> – Observe the patient position (see Chapter 9.3 on page 9-2). – Elevate the patient's upper torso and continue the treatment. – Press the  button to terminate the current treatment phase and to switch to the next treatment phase.
94	Solution temperature too high	<ul style="list-style-type: none"> – Check patient line, catheter extension and solution bag lines for kinks and closed clamps. – The solution bags used must have a temperature of 15–35 °C. – The <i>sleep•safe harmony</i> should not be exposed to direct sunlight. – If the problem occurs repeatedly, terminate the treatment.
96	Air detected	<ul style="list-style-type: none"> – Outflow: Check the patient line for leakage. – Inflow: Check bag line for leakage. – Air detected in solution bag: Suspend the solution bag with the solution bag line pointing down and continue the treatment.
100	Pumps fail to deliver	<ul style="list-style-type: none"> – The treatment has been interrupted for more than 10 minutes. Press the "Continue" button to continue the treatment.
109	Patient safety clamp error	<ul style="list-style-type: none"> – Check that the patient line is correctly positioned in the line guide of the loading tray and replace the <i>sleep•safe</i> Set, if necessary.
134 134T	Pressure in drain line too high	<ul style="list-style-type: none"> – Check the drain line for kinks and closed clamps. – The outlet of the drain line must never be more than 2 meters below or above the <i>sleep•safe harmony</i>.



Warning

- If the treatment **cannot** be continued, the *sleep•safe harmony* must be switched off with the power switch.

5.5 Premature termination of a treatment after a system error / device fault



Warning

Risk of infection

- If the treatment is stopped due to an alarm (system error / device fault), follow the instructions of the attending physician.
- Use aseptic technique when disconnecting the patient connector.
- Wearing a face mask and hand disinfection is recommended.
- Observe the hygiene practices of the dialysis center and the hygiene regulations in force.

The attending physician is responsible for proper patient training.

Additional screen layout information (see Chapter 5.1 on page 5-1).

Additional screen message information (see Chapter 5.4 on page 5-6).

5.6 Emergency shutdown

If the *sleep•safe harmony* no longer responds to screen commands, an emergency shutdown can be performed.




- To do this, first flip the power switch on the rear of the *sleep-safe harmony* to off (position "0")



- Then simultaneously press the and keys for 5 seconds.

The *sleep-safe harmony* will shut down.

5.7 Power failure



This screen will be displayed if the power supply for the *sleep•safe harmony* has failed.


The time remaining until the *sleep•safe harmony* will shut down is displayed at the bottom of the screen.

Restore power.

If power is not restored within the following 10 minutes, an audible signal will be sounded.

➤ Press the  key to silence the audible signal.





If power could not be restored in time, the treatment must be stopped.

- Disconnect the patient (see Chapter 4.4.1 on page 4-28) and remove the *sleep•safe* Set (see Chapter 4.5.3 on page 4-46).

5.8 Screen failure



Warning

- If the *sleep•safe harmony* can no longer be operated via the screen, follow the instructions of the attending physician.
- Use aseptic technique when disconnecting the patient connector.
- Wearing a face mask and hand disinfection is recommended.
- Observe the hygiene practices of the dialysis center and the hygiene regulations in force.

The attending physician is responsible for proper patient training.

6 Cleaning / disinfection

6.1 Cleaning



Warning

- Do not autoclave or submerge the *sleep•safe harmony*.
- Do not use solvent-based chemical cleaning agents.
- Disconnect the *sleep•safe harmony* from the power supply before cleaning.
- Use the recommended disinfectant and cleaning agent.

The test procedure by which the efficiency of each required disinfection has been verified is available on request.

Approved disinfectants

- Fresenius ClearSurf
- Fresenius Freka-NOL

Cleaning interval

- As necessary (in case of contamination), but at least once a week

Parts to be cleaned

- Surfaces (housing, screen, loading tray)

Once the disinfectant and cleaning agent has completely evaporated, the *sleep•safe harmony* is ready for use again.

7 Functional description

7.1 Description of functional procedures

The *sleep•safe harmony* is a peritoneal dialysis device which has been designed to offer maximum safety and convenience for patients, physicians, nurses and technicians. It reflects the latest state of technology in the fields of electronics, mechanics and software.

The most important features of the *sleep•safe harmony* are:

- Active dialysate pumps
- High flow rates
- Automatic bag detection and connection
- 100 % PVC-free line sets made of Biofine®
- Inline flow heating (no need to preheat the solution bag)
- Graphic user guidance (color touchscreen)
- Extremely quiet operation
- Integrated management system for patient and treatment data using a patient card

The *sleep•safe harmony* is an automated peritoneal dialysis device designed for use at home and in hospitals. The pumping action is enabled by a hydraulically-driven membrane piston pump.

The *sleep•safe* Set is a single use item with semi-spherical chambers in which the hydraulic system either forces dialysis solution out or draws dialysis solution into the chamber through a membrane. A built-in heater warms the dialysis solution before it is infused into the patient. An automatic connecting system connects the solution bags to the line set.


Several pressure sensors monitor the pressures. The *sleep•safe harmony* is controlled and monitored by a microprocessor. The *sleep•safe harmony* is provided with independent safety systems.

Further information material is available from the local service support organization in the form of manuals, posters etc.

7.2 Therapy types

7.2.1 Standard prescription

The following schematic shows the cycles of a standard prescription.



- 1 Treatment cycle
- 2 Initial outflow
- 3 Inflow volume
- 4 Dwell duration
- 5 Outflow volume
- 6 Last inflow
- 7 Permitted residual volume
- 8 Permitted patient volume


The standard prescription consists of individual base cycles with the possibility of profiling the inflow volume, the solution and the dwell duration per cycle.

The following parameters can be set:

- Initial outflow volume
- Number of base cycles
- Inflow volume (can be profiled)
- Solution (can be profiled)
- Dwell duration (can be profiled)
- Last inflow volume
- Last inflow solution
- Prescription name

7.2.2 Adapted APD prescription

The following schematic shows the cycles of the adapted APD prescription.



1 Cycle 1–2

2 Initial outflow

3 Cycle 1–2 inflow volume

4 Cycle 1–2 dwell duration

5 Cycle 1–2 outflow volume

6 Cycle 3–5

- 7** Cycle 3-5 inflow volume
- 8** Cycle 3-5 dwell duration
- 9** Cycle 3-5 outflow volume
- 10** Last inflow
- 11** Permitted residual volume
- 12** Permitted patient volume


The adapted APD prescription consists of two cycles with a low inflow volume and short dwell duration, and three cycles with a higher inflow volume and a longer dwell duration.

The following parameters can be set:

- Initial outflow volume
- Cycle 1–2 inflow volume
- Cycle 1–2 solution
- Cycle 1–2 dwell duration
- Cycle 3–5 inflow volume
- Cycle 3–5 solution
- Cycle 3–5 dwell duration
- Last inflow volume
- Last inflow solution
- Prescription name

7.2.3 Tidal prescription

The following schematic shows the cycles of a Tidal prescription.



- 1 Base cycle
- 2 Tidal cycles
- 3 Initial outflow
- 4 Base inflow volume including first Tidal inflow volume
- 5 Tidal inflow volume
- 6 Dwell duration
- 7 Tidal outflow volume
- 8 Base outflow volume including last Tidal outflow volume
- 9 Last inflow
- 10 Permitted residual volume
- 11 Permitted patient volume


The Tidal prescription consists of individual base cycles and Tidal cycles. The base inflow solution, the Tidal inflow solution and the dwell duration can be profiled.

The following parameters can be set:

- Initial outflow volume
- Number of base cycles
- Base inflow volume
- Number of Tidal cycles
- Tidal inflow volume
- Tidal outflow volume
- Solution (can be profiled)
- Dwell duration (can be profiled)
- Last inflow volume
- Last inflow solution
- Name

7.2.4 Basic prescription

The following schematic shows the cycles of a basic prescription.



- 5 Outflow volume
- 6 Last inflow
- 7 Permitted residual volume
- 8 Permitted patient volume


The number of cycles of a base prescription depends on the maximum inflow volume, the treatment volume and the total treatment duration.

The following parameters can be set:

- Initial outflow volume
- Maximum inflow volume
- Treatment duration
- Treatment volume
- Base cycle solution
- Last inflow volume
- Last inflow solution
- Prescription name

7.2.5 PD-Plus prescription

The following schematic shows a PD-Plus prescription.



- 2** Initial outflow
- 3** PD-Plus inflow
- 4** PD-Plus pause duration
- 5** PD-Plus outflow
- 6** Normal night-time treatment
- 7** Permitted residual volume
- 8** Permitted patient volume

The PD-Plus prescription consists of the normal night-time treatment and an additional daytime outflow and inflow. Inflow volume, solution and dwell duration can be profiled.

The following parameters can be set:

- Initial outflow volume
- Number of base cycles of the nighttime treatment
- PD-Plus inflow volume
- PD-Plus solution
- PD-Plus pause duration
- Inflow volume (can be profiled)
- Solution (can be profiled)
- Dwell duration (can be profiled)
- Last inflow volume
- Last inflow solution
- Prescription name

7.3 Therapy options

The following options must be individually adapted to the respective patient and can be set by specialized medical staff in the Patient options menu. They influence the volume and time management of the treatment.

7.3.1 Volume optimization

7.3.1.1 Permitted patient volume

Due to individual patient-related circumstances, an outflow may occasionally not be fully completed. As a consequence, a certain residual volume will remain in the peritoneal cavity. In such cases, the "permitted patient volume" option will be applied.

The permitted patient volume restricts the maximum volume which is allowed to be in the patient's peritoneal cavity over the entire treatment.

The maximum inflow volume of a prescription is computed against the percentage factor of the "permitted patient volume", resulting in the permitted patient volume for the entire treatment. Assuming a maximum inflow volume of 2000 ml and a factor of 110 % for the "permitted patient volume", the permitted patient volume would be 2200 ml.

The goal of this option is to attain the inflow volume prescribed by the physician as closely as possible, to ensure the efficacy of the dialysis treatment.

7.3.1.2 Permitted residual volume

The "permitted residual volume" is closely related to the "permitted patient volume" described above.

The permitted residual volume describes the maximum residual volume which is allowed to remain in the patient's peritoneal cavity if the outflow was not fully completed, before the device is permitted to switch to the next inflow without generating an alarm.

The goal of this option is to adapt the device to the individual catheter performance of the patient.

7.3.1.3 Permitted reduction of the inflow volume

Under certain circumstances, the prescribed total solution volume may not be available for the treatment (due to the solution volume required to prime the line set, for example). In such cases, the inflow volume will be slightly reduced.

Assuming an inflow volume of 1000 ml and a factor of 10 % for the "permitted reduction of the inflow volume", the inflow volume may be reduced by up to 100 ml, as required. The volume will only be reduced as required and within this specified limit, and the reduction is performed automatically by the device. This ensures that the treatment can be terminated without an alarm even if the originally prescribed treatment volume was not reached.

7.3.2 Time optimization

7.3.2.1 Permitted reduction of the dwell duration

The goal of the "permitted reduction of the dwell duration" option is to keep to the prescribed total treatment duration as precisely as possible. This is achieved by a dynamic adaption of the dwell duration over the remaining cycles. Assuming a dwell duration of 100 minutes and a factor of 15 % for the "permitted reduction of the dwell duration", this dwell duration may be reduced by 15 minutes.

In case of frequent reductions of the dwell duration, the parameter **Catheter performance** can be used to influence the expected treatment duration. The **Catheter performance** parameter can be changed by specialized clinical staff.

7.3.2.2 Catheter performance

The catheter performance describes the patient's outflow characteristics.

It can be set by specialized clinical staff in the **Patient parameters** menu.

The **catheter performance** influences the calculated outflow duration and thus the expected total treatment duration.

In case of frequent reductions of the dwell duration, the expected treatment duration can be adapted to the actual treatment duration using the parameter **Catheter performance**. Increasing this parameter will increase the expected treatment duration.

8 Consumables, accessories, additional equipment

The dialysis solutions required for the treatment must comply with the German Medicinal Products Act (Arzneimittelgesetz, AMG).



Warning

The device has been approved for use with certain consumables and accessories (see Chapter 8 on page 8-1). Should the responsible organization wish to use other consumables and accessories than those listed in the relevant chapter, their suitability must be verified beforehand by obtaining the appropriate information from the manufacturers.

The applicable legal regulations must be complied with.

The manufacturer does not assume any responsibility or liability for personal injury or other damage and excludes any warranty for damage to the device resulting from the use of non-approved or unsuitable consumables or accessories.

On request, the local service support organization will provide information on further accessories, consumables and other additional equipment.



Note

Various consumables are required to perform a treatment. When using dialysis solutions, make sure these solutions are used in accordance with the manufacturer's specifications.



Warning

Observe the notes in the instructions for use of the consumables concerned.



Note

Consumables:

When using consumables, it is important to take note of the following symbols:

Single use item

Identified by the symbol:



Do not reuse.

Use by date

Identified by the symbol:



Use by:

8.1 Consumables

8.1.1 Dialysis solutions

For detailed information on the dialysis solutions, refer to the extensive product range of Fresenius Medical Care.

Only dialysis solutions which are displayed on the screen under **Select solution** can be used for the treatment (see Chapter 3.4.1 on page 3-11).

8.1.2 Single use items

Product	Part number	Information	
<i>sleep•safe</i> Set	501 692 1	Line set	
<i>sleep•safe</i> Set Plus	501 693 1	Line set with second patient connector	
Drainage Extension Line - 12 m Biofine®	F00001511	Drain line extension made of Biofine®, length 12 m	
Drainage Extension Line - 12 m PVC	501 461 1	Drain line extension made of PVC, length 12 m	
Safe•Lock® PD-NIGHT® Drainage Set	501 957 1	Drain line set with 2 drain bags	
<i>stay•safe</i> ® disinfection cap	284 509 1	Disinfection cap	
APD sample bag	505 508 1	Sampling bag	
<i>sleep•safe</i> Luer lock adapter	501 697 1	Only in combination with a 7.5 % polyglucose solution bag	
Dialysis solution <i>sleep•safe balance</i> 5000 ml x 2	D, CH, A, GB, NL, E, P, GR, CY, F, B, I, S, FIN, N, DK	F 000 01 672	1.5 % glucose 1.25 mmol/L Ca
	RUS, LV, RO, LT, PL, EST, SLO, H, SK, HR+BiH, TR, CZ	F 000 00 437	

8.1.3 Surface disinfection / surface cleaning

Product	Information
Fresenius Medical Care ClearSurf	Disinfectant for wipe disinfection Active substance base: cationic surfactants
Fresenius Medical Care ClearSurf Wipes	Ready to use disinfection wipes, soaked with 1 % ClearSurf Active substance base: cationic surfactants
Fresenius Medical Care Freka-NOL	Quick-acting disinfectant for wipe disinfection in combination with Freka-Wipes disposable wipes Active substance: 45 % ethanol

8.2 Accessories

Product	Part number	Information
<i>sleep•safe harmony PatientCard^{Plus}</i>	M45 129 1	

8.3 Additional equipment

Product	Part number	Information
<i>sleep•safe transport case</i>	M39 023 1	

8.4 *sleep•safe harmony*

Product	Part number	Information
<i>sleep•safe harmony</i>	M20 600 1	

9 Installation

9.1 Connection requirements

9.1.1 Environment

The manufacturer has specified the device for operation in rooms that are suitable for peritoneal dialysis located in professional health care institutions, or for the home healthcare environment.

9.1.2 Power supply (electrical power network)

9.1.2.1 General requirements

The national standards and regulations must be observed when connecting the device to the power supply.

Power cable

If the power cable needs to be replaced, use only the original power cable listed in the spare parts catalog.

Rechargeable battery

Maintenance of the built-in rechargeable battery

Connect the *sleep•safe harmony* to the power supply every 6 months, use the power switch to switch the device on and leave on for 24 hours in Ready mode.

9.2 Installation requirements



Note

Variations in temperature during transport may cause water condensation on electrical parts. In the event of major variations in temperature, allow 3 hours for the device to adjust to the ambient temperature with the power switched off.

Operating conditions

Operating temperature range:

15 °C to 35 °C

Atmospheric pressure:

750 hPa (approx. altitude 2500 m) up to 1060 hPa

Relative humidity:

30 % to 75 %, temporarily 95 % (non-condensing)

Electromagnetic radiation

Do not use devices that emit any form of electromagnetic radiation (e.g., walkie-talkies, mobile phones, radio transmitters) in the vicinity of a *sleep•safe harmony* in operation. This may cause a malfunction of the *sleep•safe harmony* (see Chapter 12.6 on page 12-4).

9.3 Setting up the *sleep•safe harmony*




Note

➤ Remove the *sleep•safe harmony* from its transport case and place it on a flat, stable surface. Make sure the *sleep•safe harmony* is horizontal. The maximum incline allowed is less than 5°.

Note

➤ During treatment, the *sleep•safe harmony* should be positioned at the same level as the patient, but no more than 1 meter below this.

**Note**

The outlet of the drain line must never be more than 2 meters below or above the *sleep•safe harmony*.

9.4 Installation after shipment / transport outside of buildings

For storage, shipment and transporting the device outside of buildings, use the transport case provided (see Chapter 8.3 on page 8-4).

Follow the instructions below if the *sleep•safe harmony* is transported in its transport case, e.g., by car or railway, to a different location for operation.

9.4.1 Visual check after transport

Do not use the device if any of the following defects are detected:

- Mechanical damage to the housing
- Mechanical damage to the screen
- Mechanical damage to the card slot
- Defective power cable

➤ The device must be taken out of service.

➤ Inform the responsible organization or service support.

9.4.2 System check



Note

Variations in temperature during transport may cause water condensation on electrical parts. In the event of major variations in temperature, allow 3 hours for the device to adjust to the ambient temperature with the power switched off.

9.4.3 Extended system test



- Connect the *sleep-safe harmony* to the power supply.
- Switch the *sleep-safe harmony* on with the power switch.
- Press the  key as soon as it lights up.

After the power has been switched on, the device will first perform a functional test.




The *sleep-safe harmony* is ready for operation as soon as the screen shown on the left is displayed.

The following information will be displayed:

- The current date
- The software version
- The current time

➤ Press the nose of the smiley to move to the next operating step.

➤ If the *sleep-safe* Set is still in the *sleep-safe harmony* after a treatment has been prematurely terminated, press the button (1) on the right at the bottom of the screen to open the loading tray.




➤ Press the button to start the treatment with a patient card.

or

➤ Press the button only if you wish to start treatment without a patient card.

➤ The button can be used to replay the screen animation.



If the treatment was started without a patient card, the patient name must be confirmed.


- Press the  button only if you need to personalize the *sleep•safe harmony*. This will create a patient card with standard values.
- Press the  button to confirm the patient name.



The **Prescription** tab shows the following information:

- Patient name
- Prescription
- Maximum inflow volume
- Total treatment volume
- Expected treatment duration (shown here without time optimization)

- Press the  button to move to the next operating step.



The *sleep-safe harmony* will perform internal tests.


The loading tray opens.
During the test an audible alarm can be heard and the status indicator is red.

The next screen message will be displayed automatically.



Warning

- Do not use the *sleep-safe harmony* if the audible alarm is not heard or the status indicator is not lit during the internal test.
- Call service support.



➤ Do not place a *sleep-safe* Set in the open loading tray.


➤ Press the key.



- Press the button.
- Press the button to return to the previous operating step.



- The loading tray closes.
- The next screen message will be displayed automatically.



The *sleep-safe harmony* will shut down.

- If the "Extended system test" has successfully completed, the *sleep-safe harmony* can be used for treatment.
- If the "Extended system test" cannot be executed as described, do not perform a treatment.
Call service support.

10 Transport / storage

10.1 Transport within buildings

The *sleep•safe harmony* can be put back into operation without any additional tests if it is only briefly disconnected from the power supply and moved to a different room within the same building.

10.2 Shipment / transport outside of buildings

(see Chapter 9.4 on page 9-3)

10.3 Storage

The *sleep•safe harmony* must be stored in its transport case on a stable, horizontal surface in a well ventilated room.

10.3.1 Storage conditions

Temperature	-10 to +60 °C
Relative humidity	30 % to 75 %, temporarily 95 % (non-condensing)
Atmospheric pressure	500 (approx. 5000 m) up to 1060 hPa
Maintenance of the built-in rechargeable battery	Connect the <i>sleep•safe harmony</i> to the power supply every 6 months, use the power switch to switch the device on and leave on for 24 hours in Ready mode.

10.4 Environmental compatibility / disposal

Within the EC member states, the device must be disposed of in accordance with the "Directive on waste electrical and electronic equipment" (WEEE directive). Please also observe the applicable local regulations.

Before the device is sent off for disposal, the responsible organization must ensure that all consumables attached to the device are removed and the device is disinfected as specified by the manufacturer (see Chapter 6 on page 6-1).

Moreover, the responsible organization must ensure that the waste disposal company is informed of the following facts before the dismantling process is begun:

- There is a potential risk that the device is contaminated when it is returned. Suitable safety measures must be in place for dismantling the device, including the mandatory use of personal protective equipment.
- For information on the batteries and other materials used, consult the appropriate chapters of these instructions for use (see Chapter 12.10 on page 12-11) or (see Chapter 12.12 on page 12-14).
- Batteries must be properly disposed of in accordance with the applicable national regulations.
- The *sleep•safe harmony* contains triacetin, a hydraulic fluid. Triacetin must not be disposed of with domestic refuse and must not enter the sewage system. The disposal options depend on the local, national regulations and requirements, which must be followed.
- The *sleep•safe harmony* has an 8.4" TFT LC display.
- More information will be made available by the manufacturer to waste disposal services on request.

11 Technical safety checks/maintenance procedures

11.1 Important information on the technical safety checks/maintenance procedures

Technical safety checks (TSCs)	The first TSC is required before the end of the 24th month following initial start-up after delivery from the factory. All further TSCs are required before the end of the 24th month following the last TSC performed.
Maintenance procedures (MA)	The maintenance procedures (MA) are a recommendation of the manufacturer. The maintenance procedures help ensure trouble-free operation, and must be carried out for the first time before the end of the 24th month following initial start-up after delivery from the factory. All further maintenance procedures should be performed before the end of the 24th month following the last maintenance procedure performed.
Qualification requirements of testers	The checks must be performed by the Fresenius Medical Care service support organization or a person authorized by them.
	The checks must be performed by personnel qualified to perform them correctly, based on their education, training, knowledge and experience. Furthermore, the persons performing the checks must be permitted to do so independently and without outside interference.
Specifications	The information contained in the Specifications chapter must be observed.
Documentation	The Technical Safety Checks (TSC), maintenance procedures and detailed explanations of how to perform them are described in the Service Manual. Reports can be supplied on request.

The completion of the TSCs must be entered in the Medical Device Register.


12 Specifications

12.1 Dimensions and weight

Dimensions	Height: 30.5 cm Width: 45.5 cm Depth: 42.5 cm
Weight	Approx. 20 kg

12.2 Identification label (*sleep•safe harmony* product marking)

The identification label shown is only an example. Always go by the information shown on the identification label affixed to the *sleep•safe harmony* device itself.



- 1 Type identification, serial number
- 2 Power requirements
- 3 Symbol for the marking of electrical and electronic equipment
- 4 CE mark
- 5 Manufacturer and year of manufacture

- 6 Protection against ingress of solid objects:
protected against the ingress of foreign bodies larger than 12 mm in diameter; finger protection
Protection against ingress of fluids:
drip-proof
- 7 Degree of protection against electric shock: Type BF
- 8 Type of protection against electric shock:
Protection class II

12.3 Electrical safety

Classification according to EN 60601-1, IEC 60601-1

Type of protection against electric shock	Protection class II, symbol: 
Degree of protection against electric shock	Type BF, symbol 
Applied part	<i>sleep•safe</i> Set
Degree of protection against ingress of foreign bodies and liquids	IP21, symbol: IP21 <ul style="list-style-type: none">– protected against the ingress of foreign bodies larger than 12 mm in diameter; finger protection– drip-proof
Leakage currents	According to EN 60601-1

12.4 Electrical power supply

Line voltage	100–240 V AC, 50/60 Hz (Always go by the line voltage, frequency and current consumption information specified on the identification label affixed to the <i>sleep•safe harmony</i> device itself.)
Power supply connection	According to local regulations for electrical power supply.

Operating current	5 A (at 240 V AC) 8 A (at 100 V AC)
Power supply (internal)	+24 V $\pm 5\%$, 3 A, short-circuit-proof +12 V $\pm 5\%$, 0.5 A, short-circuit-proof +5 V $\pm 5\%$, 3 A, short-circuit-proof
Power switch	All-pole, simultaneous disconnection

12.5 Fuses

PCB/ module	Item no.	Value
LP1365	F1	T 0.315 A / 305 VAC / 250 VDC
LP1365	F2	T 4.0 A / 305 VAC / 250 VDC
LP1365	F3	T 8.0 A / 250 VAC / 250 VDC
LP1365	F4	T 0.125 A / 305 VAC / 250 VDC
LP1366	F1	T 3.0 A
LP1366	F2	T 1.5 A
LP1366	F3	T 5.0 A
LP1366	F4	T 4.0 A
LP1366	F5	T 3.0 A
LP1366	F6	T 1.0 A
LP1367	F1	F 1 A / 125 V
LP1367	F2	F 1 A / 125 V
LP1367	F3	F 1 A / 125 V
LP1367	F4	F 1 A / 125 V
LP1367	F5	F 1 A / 125 V

PCB/ module	Item no.	Value
LP1367	F6	F 1 A / 125 V
LP1367	F7	0.160 A / 125 V
LP1367	F8	T 3.0 A / 125 V

12.6 Information on electromagnetic compatibility

Medical electrical devices are subject to special protective measures with regard to electromagnetic compatibility (EMC).

Portable and mobile high-frequency RF communication equipment is a source of radiated electromagnetic interference that can affect medical electrical equipment. This can lead to device malfunctions.

For this reason, a certain minimum distance must be kept between such radiation sources and running medical electrical equipment.



Warning

Patient hazard due to device malfunctions

A device malfunction can be caused by high-frequency electromagnetic waves (radiated interference source).

➤ Observe the minimum required distances to radiated interference sources.

- **Minimum distances between radiation source and medical electrical equipment**

RF transmission technology employed	Minimum distance to medical electrical equipment, including all connecting cables	Radiated interference source (examples) (One device/interference source may use more than one RF transmission technology)
GSM (Global System for Mobile Communications)	1.0 m	Mobile, cell phone, smartphone, tablet computer
UMTS (Universal Mobile Telecommunications System)		
LTE (Long-Term Evolution)		
DECT (Digital Enhanced Cordless Telecommunications)	0.35 m	Cordless phone
WLAN (Wireless Local Area Network)	0.22 m	Laptop/notebook, desktop computer, E-book reader WLAN repeater, WLAN router, WLAN access point, WLAN print server
Bluetooth	0.07 m	Wireless keyboard, wireless mouse, wireless speakers
Radio	0.25 m	Radio remote control (except radio remote controls by FMC)

12.6.1 EMC guidance and manufacturer's declaration

Specifications refer to the requirements of IEC 60601-1-2:2007.

- **Electromagnetic emissions**

Guidance and manufacturer's declaration – electromagnetic emissions		
The <i>sleep•safe harmony</i> is intended for use in the electromagnetic environment specified below. The customer or the user of the <i>sleep•safe harmony</i> should make sure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions according to CISPR 11	Group 1	The <i>sleep•safe harmony</i> uses RF energy solely for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions according to CISPR 11	Class B	The <i>sleep•safe harmony</i> is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions according to IEC 61000-3-3	Compliant	

- **Electromagnetic immunity**

Guidance and manufacturer's declaration – electromagnetic immunity			
The <i>sleep•safe harmony</i> is intended for use in the electromagnetic environment specified below. The customer or the user of the <i>sleep•safe harmony</i> should make sure that it is used in such an environment.			
Immunity testing	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance

Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic materials, the relative humidity should be at least 30 %.
Electrical fast transient/burst according to IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines Not applicable	The power supply quality should be that of a typical commercial and/or hospital environment.
Surges according to IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	The power supply quality should be that of a typical commercial and/or hospital environment.
Voltage dips, short interruptions and voltage fluctuations on power supply lines according to IEC 61000-4-11	< 5 % U_T (> 95 % dip in U_T) for 0.5 cycles 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (> 95 % dip in U_T) for 5 s	< 5 % U_T (> 95 % dip in U_T) for 0.5 cycles 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (> 95 % dip in U_T) for 5 s	In the event of power supply interruptions, the battery of the <i>sleep•safe harmony</i> takes over the supply without delay.
Power frequency (50/60 Hz) magnetic field according to IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be within the limits characteristic of a typical commercial or hospital environment.
Note: U_T is the AC power supply voltage prior to the application of the test level			

			<p>Portable and mobile RF communications equipment should be used no closer to any part of the <i>sleep•safe harmony</i>, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.17 \sqrt{P}$ <p>for 150 kHz to < 80 MHz</p> $d = 0.35 \sqrt{P}$ <p>for 80 MHz to < 800 MHz</p> $d = 0.7 \sqrt{P}$ <p>for 800 MHz to 2.5 GHz</p> <p>where P is the output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.</p> <p> Interference may occur in the vicinity of equipment marked with the following symbol.</p>
Conducted RF according to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	
Radiated RF according to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and by reflection from structures, objects and people.

a

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be theoretically predicted with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the *sleep-safe harmony* is used exceeds the applicable RF compliance level above, the *sleep-safe harmony* should be observed to verify normal operation. If any irregularities in performance are observed, additional measures may be necessary, such as reorienting or relocating the *sleep-safe harmony*.

b

Field strengths should be less than 3 V/m over the entire frequency band from 150 kHz to 80 MHz.

- **Recommended separation distances between portable and mobile RF telecommunications equipment and the *sleep-safe harmony***

Recommended separation distances between portable and mobile RF telecommunications equipment and the <i>sleep-safe harmony</i>			
Nominal power of transmitter	Separation distance according to frequency of transmitter in m		
	150 kHz to < 80 MHz	80 MHz to < 800 MHz	800 MHz to 2.5 GHz

The *sleep-safe harmony* is intended for use in an electromagnetic environment in which RF disturbances are controlled. The customer or the user of the *sleep-safe harmony* can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the *sleep-safe harmony* as recommended below, according to the maximum output power of the communications equipment.

Nominal power of transmitter	Separation distance according to frequency of transmitter in m		
	150 kHz to < 80 MHz $d = 1.17 \sqrt{P}$	80 MHz to < 800 MHz $d = 1.17 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.34 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.2	1.2	2.3
10	3.7	3.7	7.4
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and by reflection from structures, objects and people.

12.7 Operating conditions

Operating temperature range	+15 to +35 °C
Atmospheric pressure	750 (approx. altitude 2500 m) up to 1060 hPa
Relative humidity	30 % to 75 %, temporarily 95 % (non-condensing)
Stability	maximum incline allowed: 5°

12.8 Storage conditions

Temperature	-10 to +60 °C
Relative humidity	30 % to 75 %, temporarily 95 % (non-condensing)
Atmospheric pressure	500 (approx. 5000 m) up to 1060 hPa
Maintenance of the built-in rechargeable battery	Connect the <i>sleep•safe harmony</i> to the power supply every 6 months, use the power switch to switch the device on and leave on for 24 hours in Ready mode.

12.9 External connection options

Any additional equipment connected to this device must comply with the relevant IEC or ISO standards (e.g., IEC 60950-1 for information technology equipment).

Furthermore, all device configurations must comply with the requirements for medical electrical systems (see EN 60601-1:2006 section 16 and annex I).

Connecting the device to an IT network that contains components not installed and validated by the device manufacturer can introduce unknown risks for patients, operators or third parties. These risks must be identified, analyzed, evaluated and monitored by the responsible organization. For assistance, consult IEC 80001-1:2010 and annexes H5 and H6 of EN 60601-1:2006.

Any modification to an IT network that has been installed and validated by the device manufacturer can introduce new risks and therefore require a repeat analysis.

Especially problematic activities:

- Changes to the IT network configuration
- Connecting additional components and devices to the IT network
- Removing components and devices from the IT network
- Updating or upgrading components and devices in the IT network

Note that local laws take priority over the above-mentioned normative requirements. Please address any queries to the local service support organization.

LAN Intended for future use

Alarm output Function currently not available:
For connecting an external alarm indicator (nurse call).
(Volt-free alarm output, alternating contact, maximum 24 V/24 W)
Port: 3-pin diode plug via a shielded line.

12.10 Batteries

The operator has no access to the batteries or the rechargeable battery.

Batteries ETX baseboard LP1450 / lithium CR 2032

Rechargeable battery Lead-acid battery (maintenance-free)
12 V (2 x 6 V), 1.3 Ah

12.11 Parameters

STANDARD inflow volume <i>sleep•safe set</i>		Adjustment range	25–3500 ml	
		Resolution	1 ml	
		Tolerance	±5 %	
TIDAL inflow volume <i>sleep•safe set</i>	base inflow volume	Adjustment range	50–3500 ml	
		Resolution	1 ml	
		Tolerance	±5 %	
	Tidal inflow volume	Adjustment range	25-3475 ml	
		Resolution	1 ml	
		Tolerance	±5 %	
Dwell duration		Adjustment range	25–3500 ml	
		Resolution	1 ml	
		Tolerance	±5 %	
STANDARD number of cycles		Adjustment range	1-50	
		Resolution	1	
TIDAL number of cycles	Base cycles	Adjustment range	1-5	
		Resolution	1	
	Tidal cycles	Adjustment range	2-15	
		Resolution	1	
Inflow rate <i>sleep•safe set</i>		Adjustment range	50–350 ml/min	
		Resolution	1 ml/min	
Outflow rate <i>sleep•safe set</i>		Adjustment range	50-230 ml/min	
		Resolution	1 ml/min	

Inflow time	Depending on flow rate	
Outflow time	Depending on flow rate	
Dosing tolerance	±3 % if dialysis solution is free of air	
Balancing tolerance	±1 % if dialysis solution is free of air	
Heater	Fixed patient inflow temperature	37 °C
	Tolerance	-7 to +2 °C
Outflow pressure monitored by the safety system	Max. -100 mbar at patient access with closed patient line	
Technique and sensitivity of the safety system against infusion of air	Pressure / volume method 6 % air per pump chamber	
Technique and sensitivity of the safety system against overfilling	Volume balancing 10 % above maximum permitted inflow volume	
Audible signal silencing	The audible signal can be silenced for 6 minutes.	
Permitted patient volume	Adjustment range 100–120 % of the maximum prescribed inflow volume,	Default value 110 %
Permitted residual volume	Adjustment range 10–50 % of the maximum volume in the patient	Default value 35 %
Permitted reduction of the dwell duration	Adjustment range 0–30 % of the prescribed dwell duration	Default value 15 %
Permitted reduction of the inflow volume	Adjustment range 10–30 % of the prescribed inflow volume	Default value 10 %
Catheter performance	Adjustment range 100–130 % = High 131–160 % = Normal 161–200 % = Low	Default value 130 %

Temperature monitoring	The inflow temperature is monitored by a temperature sensor. The limit is set to 40 °C. This value has a tolerance of ± 1 °C. In the event of an alarm situation, the alarm will be given after a maximum of 10 seconds.
Air detection	Continuous air detection during each cycle. Larger amounts of air in the inflow phase will be diverted to the drain. The operator will be alerted if air is drawn in directly.

12.12 Materials used

- Plastics and elastomers**

Abbreviation	Material	Used in
ABS	ABS (TSG) Rotec	Housing
ZK/CR	Cellular rubber	Sound insulation
POM	Delrin	Hydraulic pump housing
PA (PA6.6)	Polyamide	Barcode scanner holder
PE	Polyethylene	Container (hydraulic pump) Transport case, foam insert for transport case
PES	Polyethersulfone	Heater frame
EPDM	Ethylene propylene diene monomer	Seals, air cushion
PTFE	Polytetrafluoroethylene	Slide bushings
PVC	Polyvinyl chloride	Cable insulation
PP	Polypropylene	Wheels of transport case
	Kapton film	Transformer/electromagnet
	Silicone	Loading tray and machine block coating, cable insulation, seals

Abbreviation	Material	Used in
	Butyl rubber	Dome membrane
	Epoxy fiberglass	Printed circuit boards, insulating boards
	Keratherm® - Softtherm®	Heat conductor pad for ETX CPU

- **Metals**

Abbreviation	Material	Usage
	Aluminum	Loading tray Machine block, function block and PSU chassis Compressor retaining plate Telescopic handle for transport case
	Iron	Transformer/electromagnet, ferrite core
	Steel	Sockets, tubes, shielding plates
	Copper	Transformer/electromagnet, plug connectors, cables
	Tin	Plug connectors

- **Glass / ceramics**

Abbreviation	Material	Usage
	Borosilicate glass	Glass cylinder of hydraulic pump, screen plate of display panel
	Ceramics (aluminum oxide, Al ₂ O ₃)	Supports for heating conductors
	Ceramics (aluminum nitride AlN)	Mounting plate for heating element

- **Batteries**

Abbreviation	Material	Usage
	Lithium battery CR2032	LP1450
	Lead-acid battery, 6 V 1.3 Ah	Rechargeable battery

- **Auxiliary materials**

Abbreviation	Material	Usage
	Loctite 648, 243, 406, 603, 2701, 595, 5366 UHU-Plus Endfest 300	Bonding adhesive for screws/bolts, drive sprockets, holding frames, temperature sensors
	Enamel: Alexit 346-18 Hardener: Alexit 345-39 Undercoating: Alexit 343- 67 Hardener: Alexit 345-67	Coating system for base enclosure of the housing, lid, display panel, card slot
	Triacetin	Hydraulic fluid
	POLYLUB® GLY 151	Special lubricating grease for function block
	UNISILKON L 250 L	Special lubricating grease for function block
	Baysilone® MPH 1000	Oil for function block

12.13 Factory settings

Parameters	Factory setting	Selectable
Preselected inflow volume	2000 ml	25–3500 ml
Maximum inflow volume	3500 ml	25–3500 ml
Preselected inflow rate	350 ml/min	25–350 ml/min

Parameters	Factory setting	Selectable
Maximum inflow rate	350 ml/min	25–350 ml/min
Preselected outflow rate	230 ml/min	25-230 ml/min
Maximum outflow rate	230 ml/min	25-230 ml/min
Preselected dwell duration	120 min	5–300 min
Maximum dwell duration	300 min	5–300 min
Permitted patient volume	110 %	100–120 %
Permitted residual volume	35 %	10-50 %
Permitted dwell duration reduction	15 %	0-30 %
Permitted inflow volume reduction	10 %	10-30 %
Catheter performance	130 %	100-200 %
Set type	Standard	---
Additional outflow	Yes (with sound)	No Yes (without sound) Yes (with sound)
Access level	Select prescription	No changes Select prescription Edit prescription
Screen saver	5 min	---
Volume	100 %	---
Brightness	100 %	---

13 Definitions

13.1 Definitions and terms

Aseptic	Sterile
Autoclave	To sterilize by subjecting to high pressure saturated steam
CCPD	Continuous Cyclic Peritoneal Dialysis (continuous cyclic peritoneal dialysis)
Dialysate	The exchange fluid used for the treatment.
Hypokalemia	Potassium deficiency
Hypovolemia	Decreased blood volume
IPD	Intermittent Peritoneal Dialysis (intermittent peritoneal dialysis)
NIPD	Nightly Intermittent Peritoneal Dialysis (nightly intermittent peritoneal dialysis)
Pediatrics	Medical care of infants, children, and adolescents
Patient safe condition	Patient safety clamp closed. Dialysate flow from / to the patient stopped
Peritonitis	Inflammation of the peritoneum
Tidal	Cyclic inflow and outflow with permitted residual volume in patient's peritoneal cavity

13.2 Abbreviations

TSC	Technical Safety Checks
MA	Maintenance procedures

13.3 Symbols



Caution

IP21

Protection against ingress of solid objects:
protected against the ingress of foreign bodies larger than
12 mm in diameter; finger protection



Protection against ingress of fluids:
drip-proof



Degree of protection against electric shock:
Type BF

I/O

Type of protection against electric shock:
Protection class II

ON / OFF (power supply)



On / Off (standby)



10Base-T Ethernet port (LAN)




Batch code



Reference number



Date of manufacture



Permitted storage temperature range



Do not reuse



Latex-free



Use by



Sterile fluid path. Sterilized using steam or dry heat



Quantity



The CE mark documents compliance with the MDD 93/42/EEC.

(MDD: medical device directive)

Notified body:

TÜV PRODUCT SERVICE 0123



Symbol for the marking of electrical and electronic equipment



Alarm output



YYY

Manufacturer and year of manufacture



Serial number



Follow instructions for use

13.4 Certificates

The current versions of the EC certificates will be provided by your local service support organization on request.

14 Options

Chapter without content.

15 Appendix

15.1 *sleep-safe harmony* icons list

	On / Off key		External key for confirming actions
	Start icon		Open loading tray button
	Service login button		Repeat animation button
	Confirm input button		Discard input button
	Skip operating step button		Return to previous operating step button
	Repeat button		Return to higher-level screen button
	Patient options menu		Device options menu
	Move selection up button		Move selection down button
	Increase value button		Reduce value button
	Icon for Audible signal Off		Icon for safety-relevant message
	Icon for time		Icon for treatment duration
	Icon for patient volume		Icon for volume balance

	Icon for treatment progress		Help button
	Icon for total treatment volume		Icon for patient card inserted
	Icon for inflow		Icon for outflow
	Icon for dwell duration		Icon for cycle (inflow – dwell duration – outflow)
	Icon for standard prescription		Icon for adapted APD prescription
	Icon for Tidal prescription		Icon for basic prescription
	Icon for PD-Plus prescription		Icon for unknown prescription
	Icon for parameter not profiled		Icon for parameter profiled
	Icon for calcium value		Icon for solution (dialysate)
	Icon for glucose concentration		Icon for battery capacity
	Icon for manual outflow		Icon for please wait

15.2 Disconnection with the PIN-Reload device

PIN Reload should be used in combination with the "Pause with disconnection" option (see Chapter 4.7.2 on page 4-94).

- Before disconnecting from the device have the following ready:
 - Disinfection cap
 - PIN-Reload
 - Face mask
 - Hand disinfectant



- Turn the blue knob on the patient connector clockwise.
- Then firmly push the blue knob into the patient connector.
- Close the white clamp on the catheter extension.



- Insert the patient connector in the *sleep-safe harmony* clip.
- Place a new disinfection cap into the holder of the organizer.
- Place the PIN-Reload device into the other holder of the organizer.



- Put on the face mask.
- Disinfect your hands and dry them carefully.



- Unscrew and discard the closing cap of the new disinfection cap.
- Unscrew the catheter extension from the patient connector of the *sleep•safe* set.
- Screw the catheter extension with the PIN firmly onto the new disinfection cap.



- Unscrew the PIN-Reload device and screw it onto the patient connector.



- Pull the catheter extension with the disinfection cap straight (without turning it) out of the organizer.

15.3 Connection with the PIN-Reload device



- Ensure that the patient connector is secured in the organizer and the PIN-Reload device is still screwed onto the patient connector.
- Remove the catheter extension from your clothing.
- Wash and dry your hands thoroughly according to the instructions of the PD center.



- Place the catheter extension into the holder of the organizer.
- Put on the face mask.
- Disinfect your hands.




- Push the syringe plunger fully in to release the PIN into the patient connector.



- Unscrew the empty PIN-Reload device from the patient connector and discard.
- Visually check that the PIN was released correctly. If not, repeat the procedure with a new PIN-Reload device.



- Unscrew the catheter extension from the disinfection cap.
- Screw the catheter extension directly onto the patient connector.
- Open the white clamp on the catheter extension.



- Remove the patient connector from the organizer.
- Turn the blue knob counterclockwise to avoid an unintended release of the PIN.
- Continue your treatment.

15.4 Instructions on the use of “Free software”

Contents

- A. Peritoneal dialysis device – “Free software”**
- B. Note required according to German Medical Devices Act**
- C. Information and remarks on the free software contained in the *sleep•safe harmony***
 - I. Gnuplot Copyright
 - II. The FreeType Project LICENSE
 - III. Boost Software License
 - IV. Apache Licence
 - V. GNU GENERAL PUBLIC LICENSE Version 2
 - VI. GNU GENERAL PUBLIC LICENSE Version 3
 - VII. The LibPng License
 - VIII. The zlib License

A. Peritoneal dialysis device – “Free software”

In addition to other software, the peritoneal dialysis device contains what is called “free software” which is subject to license conditions deviating from those of the proprietary software protected for Fresenius Medical Care and their licensors.

Some of the license conditions pertaining to such free software provide that Fresenius Medical Care is authorized to distribute the peritoneal dialysis device only if the accompanying documentation contains special information and notes, supplies license conditions and/or provides the source code of such free software. Fresenius Medical Care meets these requirements by providing the copyright notices, remarks and license texts contained in Section C. below. Please note that, if such information is printed in two languages, the English version has priority.

However, the privileges granted by copyright according to Section C. and the license texts contained therein, which relate to such free software, do not include the right to make modifications to the peritoneal dialysis device and subsequently continue use of the device with these modifications. On the contrary, the German Medical Devices Act prohibits any further operation of the peritoneal dialysis device once the software contained therein has been modified, because any medical device may only be operated in the form certified. For this reason, Section B. contains an appropriate note. In such a case, Fresenius Medical Care will

stop any further technical support for the device involved. In addition, such modifications and/or manipulations may result in the extinction of warranty claims against Fresenius Medical Care or other vendors of the peritoneal dialysis device in the event a claim has arisen or might arise in respect thereto. Any utilization of the free software contained in the peritoneal dialysis device in a manner other than that required during proper operation of the system will be solely at your own risk.

Please also note that the powers listed in Section C. apply only to the “free software” mentioned therein. Any other software contained in the peritoneal dialysis device is protected by copyright for the benefit of Fresenius and their licensors and may be used only as intended for the peritoneal dialysis device.

Both the GNU General Public License (GPL) and the GNU Lesser General Public License (LGPL) are supplied with this device. You can also download these license conditions from the Internet.

For GPL, please refer to: <http://www.gnu.org/copyleft/gpl.html>

For LGPL, please refer to: <http://www.gnu.org/copyleft/lesser.html>

B. Note required according to German Medical Devices Act

This medical device has been certified jointly with the operating system software Linux version 2.6.27. Any modification to the software contained in this medical device, including the operating system software, may result in the medical device losing its conformity with the regulations of the German Medical Devices Act and in losing its right to bear the CE mark. Anyone operating a medical device without a valid CE mark according to the Medical Device Directive 93/42/EEC is liable to be prosecuted. According to Section 41 MPG (Medizinproduktegesetz, German Medical Devices Act), he or she may be sentenced to up to one year of imprisonment or may be fined. In addition, anyone modifying the software contained in this medical device or allowing such a modification will also be subject to product liability against third parties who might be injured.

C. Information and remarks on the free software contained in the peritoneal dialysis device

I. Gnuplot Copyright

Copyright (C) 1986 - 1993, 1998, 2004, 2007 Thomas Williams, Colin Kelley

Permission to use, copy, and distribute this software and its documentation for any purpose with or without fee is hereby granted, provided that the above copyright notice appear in all copies and that both that copyright notice and this permission notice appear in supporting documentation.

Permission to modify the software is granted, but not the right to distribute the complete modified source code. Modifications are to be distributed as patches to the released version. Permission to distribute binaries produced by compiling modified sources is granted, provided you

1. distribute the corresponding source modifications from the released version in the form of a patch file along with the binaries,
2. add special version identification to distinguish your version in addition to the base release version number,
3. provide your name and address as the primary contact for the support of your modified version, and
4. retain our contact information in regard to use of the base software.

Permission to distribute the released version of the source code along with corresponding source modifications in the form of a patch file is granted with same provisions 2 through 4 for binary distributions.

This software is provided "as is" without express or implied warranty to the extent permitted by applicable law.

AUTHORS

Original Software:

Thomas Williams, Colin Kelley.

Gnuplot 2.0 additions:

Russell Lang, Dave Kotz, John Campbell.

Gnuplot 3.0 additions:

Gershon Elber and many others.

Gnuplot 4.0 additions:

See list of contributors at head of this document.

II. The FreeType Project LICENSE

2006-Jan-27

Copyright 1996-2002, 2006 by David Turner, Robert Wilhelm, and Werner Lemberg

Introduction

The FreeType Project is distributed in several archive packages; some of them may contain, in addition to the FreeType font engine, various tools and contributions which rely on, or relate to, the FreeType Project.

This license applies to all files found in such packages, and which do not fall under their own explicit license. The license affects thus the FreeType font engine, the test programs, documentation and makefiles, at the very least.

This license was inspired by the BSD, Artistic, and IJG (Independent JPEG Group) licenses, which all encourage inclusion and use of free software in commercial and freeware products alike. As a consequence, its main points are that:

- We don't promise that this software works. However, we will be interested in any kind of bug reports. ('as is' distribution)
- You can use this software for whatever you want, in parts or full form, without having to pay us. ('royalty-free' usage)
- You may not pretend that you wrote this software. If you use it, or only parts of it, in a program, you must acknowledge somewhere in your documentation that you have used the FreeType code. ('credits')

We specifically permit and encourage the inclusion of this software, with or without modifications, in commercial products. We disclaim all warranties covering The FreeType Project and assume no liability related to The FreeType Project.

Finally, many people asked us for a preferred form for a credit/disclaimer to use in compliance with this license. We thus encourage you to use the following text:

"Portions of this software are copyright © <year> The FreeType Project (www.freetype.org). All rights reserved."

Please replace <year> with the value from the FreeType version you actually use.

Legal Terms

0. Definitions

Throughout this license, the terms 'package', 'FreeType Project', and 'FreeType archive' refer to the set of files originally distributed by the authors (David Turner, Robert Wilhelm, and Werner Lemberg) as the 'FreeType Project', be they named as alpha, beta or final release.

'You' refers to the licensee, or person using the project, where 'using' is a generic term including compiling the project's source code as well as linking it to form a 'program' or 'executable'. This program is referred to as 'a program using the FreeType engine'.

This license applies to all files distributed in the original FreeType Project, including all source code, binaries and documentation, unless otherwise stated in the file in its original, unmodified form as distributed in the original archive. If you are unsure whether or not a particular file is covered by this license, you must contact us to verify this.

The FreeType Project is copyright © 1996-2000 by David Turner, Robert Wilhelm, and Werner Lemberg. All rights reserved except as specified below.

1. No Warranty

THE FREETYPE PROJECT IS PROVIDED 'AS IS' WITHOUT WARRANTY OF ANY KIND, EITHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL ANY OF THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY DAMAGES CAUSED BY THE USE OR THE INABILITY TO USE, OF THE FREETYPE PROJECT.

2. Redistribution

This license grants a worldwide, royalty-free, perpetual and irrevocable right and license to use, execute, perform, compile, display, copy, create derivative works of, distribute and sublicense the FreeType Project (in both source and object code forms) and derivative works thereof for any purpose; and to authorize others to exercise some or all of the rights granted herein, subject to the following conditions:

- Redistribution of source code must retain this license file ('FTL.TXT') unaltered; any additions, deletions or changes to the original files must be clearly indicated in accompanying documentation. The copyright notices of the unaltered.

III. Boost Software License

Boost Software License - Version 1.0 - August 7th, 2003

Permission is hereby granted, free of charge, to any person or organization obtaining a copy of the software and accompanying documentation covered by this license (the "Software") to use, reproduce, display, distribute, execute, and transmit the Software, and to prepare derivative works of the Software, and to permit third-parties to whom the Software is furnished to do so, all subject to the following:

The copyright notices in the Software and this entire statement, including the above license grant, this restriction and the following disclaimer, must be included in all copies of the Software, in whole or in part, and all derivative works of the Software, unless such copies or derivative works are solely in the form of machine-executable object code generated by a source language processor.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, TITLE AND NON-INFRINGEMENT. IN NO EVENT SHALL THE COPYRIGHT HOLDERS OR ANYONE DISTRIBUTING THE SOFTWARE BE LIABLE FOR ANY DAMAGES OR OTHER LIABILITY, WHETHER IN CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

IV. Apache Licence

Apache License
Version 2.0, January 2004
<http://www.apache.org/licenses/>

TERMS AND CONDITIONS FOR USE, REPRODUCTION, AND DISTRIBUTION

1. Definitions

"License" shall mean the terms and conditions for use, reproduction, and distribution as defined by Sections 1 through 9 of this document.

"Licensor" shall mean the copyright owner or entity authorized by the copyright owner that is granting the License.

"Legal Entity" shall mean the union of the acting entity and all other entities that control, are controlled by, or are under common control with that entity. For the purposes of this definition, "control" means (i) the power, direct or indirect, to cause the direction or management of such entity, whether by contract or otherwise, or (ii) ownership of fifty percent (50 %) or more of the outstanding shares, or (iii) beneficial ownership of such entity.

"You" (or "Your") shall mean an individual or Legal Entity exercising permissions granted by this License.

"Source" form shall mean the preferred form for making modifications, including but not limited to software source code, documentation source, and configuration files.

"Object" form shall mean any form resulting from mechanical transformation or translation of a Source form, including but not limited to compiled object code, generated documentation, and conversions to other media types.

"Work" shall mean the work of authorship, whether in Source or Object form, made available under the License, as indicated by a copyright notice that is included in or attached to the work (an example is provided in the Appendix below).

"Derivative Works" shall mean any work, whether in Source or Object form, that is based on (or derived from) the Work and for which the editorial revisions, annotations, elaborations, or other modifications represent, as a whole, an original work of authorship. For the purposes of this License, Derivative Works shall not include works that remain separable from, or merely link (or bind by name) to the interfaces of, the Work and Derivative Works thereof.

"Contribution" shall mean any work of authorship, including the original version of the Work and any modifications or additions to that Work or Derivative Works thereof, that is intentionally submitted to Licensor for inclusion in the Work by the copyright owner or by an individual or Legal Entity authorized to submit on behalf of the copyright owner. For the purposes of this definition, "submitted" means any form of electronic, verbal, or written communication sent to the Licensor or its representatives, including but not limited to communication on electronic mailing lists, source code control systems, and issue tracking systems that are managed by, or on behalf of, the Licensor for the purpose of discussing and improving the Work, but excluding communication that is conspicuously marked or otherwise designated in writing by the copyright owner as "Not a Contribution."

"Contributor" shall mean Licensor and any individual or Legal Entity on behalf of whom a Contribution has been received by Licensor and subsequently incorporated within the Work.

2. Grant of Copyright License.

Subject to the terms and conditions of this License, each Contributor hereby grants to You a perpetual, worldwide, non-exclusive, no-charge, royalty-free, irrevocable copyright license to reproduce, prepare Derivative Works of, publicly display, publicly perform, sublicense, and distribute the Work and such Derivative Works in Source or Object form.

3. Grant of Patent License.

Subject to the terms and conditions of this License, each Contributor hereby grants to You a perpetual, worldwide, non-exclusive, no-charge, royalty-free, irrevocable (except as stated in this section) patent license to make, have made, use, offer to sell, sell, import, and otherwise transfer the Work, where such license applies only to those patent claims licensable by such Contributor that are necessarily infringed by their Contribution(s) alone or by combination of their Contribution(s) with the Work to which such Contribution(s) was submitted. If You institute patent litigation against any entity (including a cross-claim or counterclaim in a lawsuit) alleging that the Work or a Contribution incorporated within the Work constitutes direct or contributory patent infringement, then any patent licenses granted to You under this License for that Work shall terminate as of the date such litigation is filed.

4. Redistribution.

You may reproduce and distribute copies of the Work or Derivative Works thereof in any medium, with or without modifications, and in Source or Object form, provided that You meet the following conditions:

- a) You must give any other recipients of the Work or Derivative Works a copy of this License; and
- b) You must cause any modified files to carry prominent notices stating that You changed the files; and

- c) You must retain, in the Source form of any Derivative Works that You distribute, all copyright, patent, trademark, and attribution notices from the Source form of the Work, excluding those notices that do not pertain to any part of the Derivative Works; and
- d) If the Work includes a "NOTICE" text file as part of its distribution, then any Derivative Works that You distribute must include a readable copy of the attribution notices contained within such NOTICE file, excluding those notices that do not pertain to any part of the Derivative Works, in at least one of the following places: within a NOTICE text file distributed as part of the Derivative Works; within the Source form or documentation, if provided along with the Derivative Works; or, within a display generated by the Derivative Works, if and wherever such third-party notices normally appear. The contents of the NOTICE file are for informational purposes only and do not modify the License. You may add Your own attribution notices within Derivative Works that You distribute, alongside or as an addendum to the NOTICE text from the Work, provided that such additional attribution notices cannot be construed as modifying the License.

You may add Your own copyright statement to Your modifications and may provide additional or different license terms and conditions for use, reproduction, or distribution of Your modifications, or for any such Derivative Works as a whole, provided Your use, reproduction, and distribution of the Work otherwise complies with the conditions stated in this License.

5. Submission of Contributions.

Unless You explicitly state otherwise, any Contribution intentionally submitted for inclusion in the Work by You to the Licensor shall be under the terms and conditions of this License, without any additional terms or conditions. Notwithstanding the above, nothing herein shall supersede or modify the terms of any separate license agreement you may have executed with Licensor regarding such Contributions.

6. Trademarks.

This License does not grant permission to use the trade names, trademarks, service marks, or product names of the Licensor, except as required for reasonable and customary use in describing the origin of the Work and reproducing the content of the NOTICE file.

7. Disclaimer of Warranty.

Unless required by applicable law or agreed to in writing, Licensor provides the Work (and each Contributor provides its Contributions) on an "AS IS" BASIS, WITHOUT WARRANTIES OR CONDITIONS OF ANY KIND, either express or implied, including, without limitation, any warranties or conditions of TITLE, NON-INFRINGEMENT, MERCHANTABILITY, or FITNESS FOR A PARTICULAR PURPOSE. You are solely responsible for determining the appropriateness of using or redistributing the Work and assume any risks associated with Your exercise of permissions under this License.

8. Limitation of Liability.

In no event and under no legal theory, whether in tort (including negligence), contract, or otherwise, unless required by applicable law (such as deliberate and grossly negligent acts) or agreed to in writing, shall any Contributor be liable to You for damages, including any direct, indirect, special, incidental, or consequential damages of any character arising as a result of this License or out of the use or inability to use the Work (including but not limited to damages for loss of goodwill, work stoppage, computer failure or malfunction, or any and all other commercial damages or losses), even if such Contributor has been advised of the possibility of such damages.

9. Accepting Warranty or Additional Liability.

While redistributing the Work or Derivative Works thereof, You may choose to offer, and charge a fee for, acceptance of support, warranty, indemnity, or other liability obligations and/or rights consistent with this License. However, in accepting such obligations, You may act only on Your own behalf and on Your sole responsibility, not on behalf of any other Contributor, and only if You agree to indemnify, defend, and hold each Contributor harmless for any liability incurred by, or claims asserted against, such Contributor by reason of your accepting any such warranty or additional liability.

END OF TERMS AND CONDITIONS

APPENDIX: How to apply the Apache License to your work.

To apply the Apache License to your work, attach the following boilerplate notice, with the fields enclosed by brackets "[]" replaced with your own identifying information. (Don't include the brackets!) The text should be enclosed in the appropriate comment syntax for the file format. We also recommend that a file or class name and description of purpose be included on the same "printed page" as the copyright notice for easier identification within third-party archives.

Copyright [yyyy] [name of copyright owner]

Licensed under the Apache License, Version 2.0 (the "License"); you may not use this file except in compliance with the License. You may obtain a copy of the License at
<http://www.apache.org/licenses/LICENSE-2.0>

Unless required by applicable law or agreed to in writing, software distributed under the License is distributed on an "AS IS" BASIS, WITHOUT WARRANTIES OR CONDITIONS OF ANY KIND, either express or implied. See the License for the specific language governing permissions and limitations under the License.

V. GNU GENERAL PUBLIC LICENSE Version 2

GNU GENERAL PUBLIC LICENSE

Version 2, June 1991

Copyright (C) 1989, 1991 Free Software Foundation, Inc., 51 Franklin Street, Fifth Floor, Boston, MA 02110-1301 USA Everyone is permitted to copy and distribute verbatim copies of this license document, but changing it is not allowed.

Preamble

The licenses for most software are designed to take away your freedom to share and change it. By contrast, the GNU General Public License is intended to guarantee your freedom to share and change free software to make sure the software is free for all its users. This General Public License applies to most of the Free Software Foundation's software and to any other program whose authors commit to using it. (Some other Free Software Foundation software is covered by the GNU Lesser General Public License instead.) You can apply it to your programs, too.

When we speak of free software, we are referring to freedom, not price. Our General Public Licenses are designed to make sure that you have the freedom to distribute copies of free software (and charge for this service if you wish), that you receive source code or can get it if you want it, that you can change the software or use pieces of it in new free programs; and that you know you can do these things.

To protect your rights, we need to make restrictions that forbid anyone to deny you these rights or to ask you to surrender the rights. These restrictions translate to certain responsibilities for you if you distribute copies of the software, or if you modify it.

For example, if you distribute copies of such a program, whether gratis or for a fee, you must give the recipients all the rights that you have. You must make sure that they, too, receive or can get the source code. And you must show them these terms so they know their rights.

We protect your rights with two steps: (1) copyright the software, and (2) offer you this license which gives you legal permission to copy, distribute and/or modify the software.

Also, for each author's protection and ours, we want to make certain that everyone understands that there is no warranty for this free software. If the software is modified by someone else and passed on, we want its recipients to know that what they have is not the original, so that any problems introduced by others will not reflect on the original authors' reputations.

Finally, any free program is threatened constantly by software patents. We wish to avoid the danger that redistributors of a free program will individually obtain patent licenses, in effect making the program proprietary. To prevent this, we have made it clear that any patent must be licensed for everyone's free use or not licensed at all.

The precise terms and conditions for copying, distribution and modification follow.

GNU GENERAL PUBLIC LICENSE TERMS AND CONDITIONS FOR COPYING, DISTRIBUTION AND MODIFICATION

0. This License applies to any program or other work which contains a notice placed by the copyright holder saying it may be distributed under the terms of this General Public License. The "Program", below, refers to any such program or work, and a "work based on the Program" means either the Program or any derivative work under copyright law: that is to say, a work containing the Program or a portion of it, either verbatim or with modifications and/or translated into another language. (Hereinafter, translation is included without limitation in the term "modification".) Each licensee is addressed as "you".

Activities other than copying, distribution and modification are not covered by this License; they are outside its scope. The act of running the Program is not restricted, and the output from the Program is covered only if its contents constitute a work based on the Program (independent of having been made by running the Program). Whether that is true depends on what the Program does.

1. You may copy and distribute verbatim copies of the Program's source code as you receive it, in any medium, provided that you conspicuously and appropriately publish on each copy an appropriate copyright notice and disclaimer of warranty; keep intact all the notices that refer to this License and to the absence of any warranty; and give any other recipients of the Program a copy of this License along with the Program.

You may charge a fee for the physical act of transferring a copy, and you may at your option offer warranty protection in exchange for a fee.

2. You may modify your copy or copies of the Program or any portion of it, thus forming a work based on the Program, and copy and distribute such modifications or work under the terms of Section 1 above, provided that you also meet all of these conditions:

a) You must cause the modified files to carry prominent notices stating that you changed the files and the date of any change.

b) You must cause any work that you distribute or publish, that in whole or in part contains or is derived from the Program or any part thereof, to be licensed as a whole at no charge to all third parties under the terms of this License.

c) If the modified program normally reads commands interactively when run, you must cause it, when started running for such interactive use in the most ordinary way, to print or display an announcement including an appropriate copyright notice and a notice that there is no warranty (or else, saying that you provide a warranty) and that users may redistribute the program under these conditions, and telling the user how to view a copy of this License. (Exception: if the Program itself is interactive but does not normally print such an announcement, your work based on the Program is not required to print an announcement.)

These requirements apply to the modified work as a whole. If identifiable sections of that work are not derived from the Program, and can be reasonably considered independent and separate works in themselves, then this License, and its terms, do not apply to those sections when you distribute them as separate works. But when you distribute the same sections as part of a whole which is a work based on the Program, the distribution of the whole must be on the terms of this License, whose permissions for other licensees extend to the entire whole, and thus to each and every part regardless of who wrote it.

Thus, it is not the intent of this section to claim rights or contest your rights to work written entirely by you; rather, the intent is to exercise the right to control the distribution of derivative or collective works based on the Program.

In addition, mere aggregation of another work not based on the Program with the Program (or with a work based on the Program) on a volume of a storage or distribution medium does not bring the other work under the scope of this License.

3. You may copy and distribute the Program (or a work based on it, under Section 2) in object code or executable form under the terms of Sections 1 and 2 above provided that you also do one of the following:

- a) Accompany it with the complete corresponding machine-readable source code, which must be distributed under the terms of Sections 1 and 2 above on a medium customarily used for software interchange; or,
- b) Accompany it with a written offer, valid for at least three years, to give any third party, for a charge no more than your cost of physically performing source distribution, a complete machine-readable copy of the corresponding source code, to be distributed under the terms of Sections 1 and 2 above on a medium customarily used for software interchange; or,
- c) Accompany it with the information you received as to the offer to distribute corresponding source code. (This alternative is allowed only for noncommercial distribution and only if you received the program in object code or executable form with such an offer, in accord with Subsection b above.)

The source code for a work means the preferred form of the work for making modifications to it. For an executable work, complete source code means all the source code for all modules it contains, plus any associated interface definition files, plus the scripts used to control compilation and installation of the executable. However, as a special exception, the source code distributed need not include anything that is normally distributed (in either source or binary form) with the major components (compiler, kernel, and so on) of the operating system on which the executable runs, unless that component itself accompanies the executable.

If distribution of executable or object code is made by offering access to copy from a designated place, then offering equivalent access to copy the source code from the same place counts as distribution of the source code, even though third parties are not compelled to copy the source along with the object code.

4. You may not copy, modify, sublicense, or distribute the Program except as expressly provided under this License. Any attempt otherwise to copy, modify, sublicense or distribute the Program is void, and will automatically terminate your rights under this License. However, parties who have received copies, or rights, from you under this License will not have their licenses terminated so long as such parties remain in full compliance.

5. You are not required to accept this License, since you have not signed it. However, nothing else grants you permission to modify or distribute the Program or its derivative works. These actions are prohibited by law if you do not accept this License. Therefore, by modifying or distributing the Program (or any work based on the Program), you indicate your acceptance of this License to do so, and all its terms and conditions for copying, distributing or modifying the Program or works based on it.

6. Each time you redistribute the Program (or any work based on the Program), the recipient automatically receives a license from the original licensor to copy, distribute or modify the Program subject to these terms and conditions. You may not impose any further restrictions on the recipients' exercise of the rights granted herein. You are not responsible for enforcing compliance by third parties to this License.

7. If, as a consequence of a court judgment or allegation of patent infringement or for any other reason (not limited to patent issues), conditions are imposed on you (whether by court order, agreement or otherwise) that contradict the conditions of this License, they do not excuse you from the conditions of this License. If you cannot distribute so as to satisfy simultaneously your obligations under this License and any other pertinent obligations, then as a consequence you may not distribute the Program at all. For example, if a patent license would not permit royalty-free redistribution of the Program by all those who receive copies directly or indirectly through you, then the only way you could satisfy both it and this License would be to refrain entirely from distribution of the Program.

If any portion of this section is held invalid or unenforceable under any particular circumstance, the balance of the section is intended to apply and the section as a whole is intended to apply in other circumstances.

It is not the purpose of this section to induce you to infringe any patents or other property right claims or to contest validity of any such claims; this section has the sole purpose of protecting the integrity of the free software distribution system, which is implemented by public license practices. Many people have made generous contributions to the wide range of software distributed through that system in reliance on consistent application of that system; it is up to the author/donor to decide if he or she is willing to distribute software through any other system and a licensee cannot impose that choice.

This section is intended to make thoroughly clear what is believed to be a consequence of the rest of this License.

8. If the distribution and/or use of the Program is restricted in certain countries either by patents or by copyrighted interfaces, the original copyright holder who places the Program under this License may add an explicit geographical distribution limitation excluding those countries, so that distribution is permitted only in or among countries not thus excluded. In such case, this License incorporates the limitation as if written in the body of this License.

9. The Free Software Foundation may publish revised and/or new versions of the General Public License from time to time. Such new versions will be similar in spirit to the present version, but may differ in detail to address new problems or concerns.

Each version is given a distinguishing version number. If the Program specifies a version number of this License which applies to it and "any later version", you have the option of following the terms and conditions either of that version or of any later version published by the Free Software Foundation. If the Program does not specify a version number of this License, you may choose any version ever published by the Free Software Foundation.

10. If you wish to incorporate parts of the Program into other free programs whose distribution conditions are different, write to the author to ask for permission. For software which is copyrighted by the Free Software Foundation, write to the Free Software Foundation; we sometimes make exceptions for this. Our decision will be guided by the two goals of preserving the free status of all derivatives of our free software and of promoting the sharing and reuse of software generally.

NO WARRANTY

11. BECAUSE THE PROGRAM IS LICENSED FREE OF CHARGE, THERE IS NO WARRANTY FOR THE PROGRAM, TO THE EXTENT PERMITTED BY APPLICABLE LAW. EXCEPT WHEN OTHERWISE STATED IN WRITING THE COPYRIGHT HOLDERS AND/OR OTHER PARTIES PROVIDE THE PROGRAM "AS IS" WITHOUT WARRANTY OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE ENTIRE RISK AS TO THE QUALITY AND PERFORMANCE OF THE PROGRAM IS WITH YOU. SHOULD THE PROGRAM PROVE DEFECTIVE, YOU ASSUME THE COST OF ALL NECESSARY SERVICING, REPAIR OR CORRECTION.

12. IN NO EVENT UNLESS REQUIRED BY APPLICABLE LAW OR AGREED TO IN WRITING WILL ANY COPYRIGHT HOLDER, OR ANY OTHER PARTY WHO MAY MODIFY AND/OR REDISTRIBUTE THE PROGRAM AS PERMITTED ABOVE, BE LIABLE TO YOU FOR DAMAGES, INCLUDING ANY GENERAL, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF THE USE OR INABILITY TO USE THE PROGRAM (INCLUDING BUT NOT LIMITED TO LOSS OF DATA OR DATA BEING RENDERED INACCURATE OR LOSSES SUSTAINED BY YOU OR THIRD PARTIES OR A FAILURE OF THE PROGRAM TO OPERATE WITH ANY OTHER PROGRAMS), EVEN IF SUCH HOLDER OR OTHER PARTY HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

END OF TERMS AND CONDITIONS

How to Apply These Terms to Your New Programs

If you develop a new program, and you want it to be of the greatest possible use to the public, the best way to achieve this is to make it free software which everyone can redistribute and change under these terms.

To do so, attach the following notices to the program. It is safest to attach them to the start of each source file to most effectively convey the exclusion of warranty; and each file should have at least the "copyright" line and a pointer to where the full notice is found.

<one line to give the program's name and a brief idea of what it does.> Copyright (C) <year>
<name of author>

This program is free software; you can redistribute it and/or modify it under the terms of the GNU General Public License as published by the Free Software Foundation; either version 2 of the License, or (at your option) any later version.

This program is distributed in the hope that it will be useful, but WITHOUT ANY WARRANTY; without even the implied warranty of MERCHANTABILITY or FITNESS FOR A PARTICULAR PURPOSE. See the GNU General Public License for more details.

You should have received a copy of the GNU General Public License along with this program; if not, write to the Free Software Foundation, Inc.,
51 Franklin Street, Fifth Floor, Boston, MA 02110-1301 USA.

Also add information on how to contact you by electronic and paper mail.

If the program is interactive, make it output a short notice like this when it starts in an interactive mode:

Gnomovision version 69, Copyright (C) year name of author
Gnomovision comes with ABSOLUTELY NO WARRANTY; for details type `show w'.
This is free software, and you are welcome to redistribute it under certain conditions; type `show c' for details.

The hypothetical commands `show w' and `show c' should show the appropriate parts of the General Public License. Of course, the commands you use may be called something other than `show w' and `show c'; they could even be mouse-clicks or menu items--whatever suits your program.

You should also get your employer (if you work as a programmer) or your school, if any, to sign a "copyright disclaimer" for the program, if necessary. Here is a sample; alter the names:

Yoyodyne, Inc., hereby disclaims all copyright interest in the program
'Gnomovision' (which makes passes at compilers) written by James Hacker.

<signature of Ty Coon>, 1 April 1989
Ty Coon, President of Vice

This General Public License does not permit incorporating your program into proprietary programs. If your program is a subroutine library, you may consider it more useful to permit linking proprietary applications with the library. If this is what you want to do, use the GNU Lesser General Public License instead of this License.

VI. GNU GENERAL PUBLIC LICENSE Version 3

GNU GENERAL PUBLIC LICENSE
Version 3, 29 June 2007

Copyright (C) 2007 Free Software Foundation, Inc. <<http://fsf.org/>> Everyone is permitted to copy and distribute verbatim copies of this license document, but changing it is not allowed.

Preamble

The GNU General Public License is a free, copyleft license for software and other kinds of works.

The licenses for most software and other practical works are designed to take away your freedom to share and change the works. By contrast, the GNU General Public License is intended to guarantee your freedom to share and change all versions of a program--to make sure it remains free software for all its users. We, the Free Software Foundation, use the GNU General Public License for most of our software; it applies also to any other work released this way by its authors. You can apply it to your programs, too.

When we speak of free software, we are referring to freedom, not price. Our General Public Licenses are designed to make sure that you have the freedom to distribute copies of free software (and charge for them if you wish), that you receive source code or can get it if you want it, that you can change the software or use pieces of it in new free programs, and that you know you can do these things.

To protect your rights, we need to prevent others from denying you these rights or asking you to surrender the rights. Therefore, you have certain responsibilities if you distribute copies of the software, or if you modify it: responsibilities to respect the freedom of others.

For example, if you distribute copies of such a program, whether gratis or for a fee, you must pass on to the recipients the same freedoms that you received. You must make sure that they, too, receive or can get the source code. And you must show them these terms so they know their rights.

Developers that use the GNU GPL protect your rights with two steps:

(1) assert copyright on the software, and (2) offer you this License giving you legal permission to copy, distribute and/or modify it.

For the developers' and authors' protection, the GPL clearly explains that there is no warranty for this free software. For both users' and authors' sake, the GPL requires that modified versions be marked as changed, so that their problems will not be attributed erroneously to authors of previous versions.

Some devices are designed to deny users access to install or run modified versions of the software inside them, although the manufacturer can do so. This is fundamentally incompatible with the aim of protecting users' freedom to change the software. The systematic pattern of such abuse occurs in the area of products for individuals to use, which is precisely where it is most unacceptable. Therefore, we have designed this version of the GPL to prohibit the practice for those products. If such problems arise substantially in other domains, we stand ready to extend this provision to those domains in future versions of the GPL, as needed to protect the freedom of users.

Finally, every program is threatened constantly by software patents. States should not allow patents to restrict development and use of software on general-purpose computers, but in those that do, we wish to avoid the special danger that patents applied to a free program could make it effectively proprietary. To prevent this, the GPL assures that patents cannot be used to render the program non-free.

The precise terms and conditions for copying, distribution and modification follow.

TERMS AND CONDITIONS

0. Definitions.

"This License" refers to version 3 of the GNU General Public License.

"Copyright" also means copyright-like laws that apply to other kinds of works, such as semiconductor masks.

"The Program" refers to any copyrightable work licensed under this License. Each licensee is addressed as "you". "Licensees" and "recipients" may be individuals or organizations.

To "modify" a work means to copy from or adapt all or part of the work in a fashion requiring copyright permission, other than the making of an exact copy. The resulting work is called a "modified version" of the earlier work or a work "based on" the earlier work.

A "covered work" means either the unmodified Program or a work based on the Program.

To "propagate" a work means to do anything with it that, without permission, would make you directly or secondarily liable for infringement under applicable copyright law, except executing it on a computer or modifying a private copy. Propagation includes copying, distribution (with or without modification), making available to the public, and in some countries other activities as well.

To "convey" a work means any kind of propagation that enables other parties to make or receive copies. Mere interaction with a user through a computer network, with no transfer of a copy, is not conveying.

An interactive user interface displays "Appropriate Legal Notices" to the extent that it includes a convenient and prominently visible feature that (1) displays an appropriate copyright notice, and (2) tells the user that there is no warranty for the work (except to the extent that warranties are provided), that licensees may convey the work under this License, and how to view a copy of this License. If the interface presents a list of user commands or options, such as a menu, a prominent item in the list meets this criterion.

1. Source Code.

The "source code" for a work means the preferred form of the work for making modifications to it. "Object code" means any non-source form of a work.

A "Standard Interface" means an interface that either is an official standard defined by a recognized standards body, or, in the case of interfaces specified for a particular programming language, one that is widely used among developers working in that language.

The "System Libraries" of an executable work include anything, other than the work as a whole, that (a) is included in the normal form of packaging a Major Component, but which is not part of that Major Component, and (b) serves only to enable use of the work with that Major Component, or to implement a Standard Interface for which an implementation is available to the public in source code form. A "Major Component", in this context, means a major essential component (kernel, window system, and so on) of the specific operating system (if any) on which the executable work runs, or a compiler used to produce the work, or an object code interpreter used to run it.

The "Corresponding Source" for a work in object code form means all the source code needed to generate, install, and (for an executable work) run the object code and to modify the work, including scripts to control those activities. However, it does not include the work's System Libraries, or general-purpose tools or generally available free programs which are used unmodified in performing those activities but which are not part of the work. For example, Corresponding Source includes interface definition files associated with source files for the work, and the source code for shared libraries and dynamically linked subprograms that the work is specifically designed to require, such as by intimate data communication or control flow between those subprograms and other parts of the work.

The Corresponding Source need not include anything that users can regenerate automatically from other parts of the Corresponding Source.

The Corresponding Source for a work in source code form is that same work.

2. Basic Permissions.

All rights granted under this License are granted for the term of copyright on the Program, and are irrevocable provided the stated conditions are met. This License explicitly affirms your unlimited permission to run the unmodified Program. The output from running a covered work is covered by this License only if the output, given its content, constitutes a covered work. This License acknowledges your rights of fair use or other equivalent, as provided by copyright law.

You may make, run and propagate covered works that you do not convey, without conditions so long as your license otherwise remains in force. You may convey covered works to others for the sole purpose of having them make modifications exclusively for you, or provide you with facilities for running those works, provided that you comply with the terms of this License in conveying all material for which you do not control copyright. Those thus making or running the covered works for you must do so exclusively on your behalf, under your direction and control, on terms that prohibit them from making any copies of your copyrighted material outside their relationship with you.

Conveying under any other circumstances is permitted solely under the conditions stated below. Sublicensing is not allowed; section 10 makes it unnecessary.

3. Protecting Users' Legal Rights From Anti-Circumvention Law.

No covered work shall be deemed part of an effective technological measure under any applicable law fulfilling obligations under article 11 of the WIPO copyright treaty adopted on 20 December 1996, or similar laws prohibiting or restricting circumvention of such measures.

When you convey a covered work, you waive any legal power to forbid circumvention of technological measures to the extent such circumvention is effected by exercising rights under this License with respect to the covered work, and you disclaim any intention to limit operation or modification of the work as a means of enforcing, against the work's users, your or third parties' legal rights to forbid circumvention of technological measures.

4. Conveying Verbatim Copies.

You may convey verbatim copies of the Program's source code as you receive it, in any medium, provided that you conspicuously and appropriately publish on each copy an appropriate copyright notice; keep intact all notices stating that this License and any non-permissive terms added in accord with section 7 apply to the code; keep intact all notices of the absence of any warranty; and give all recipients a copy of this License along with the Program.

You may charge any price or no price for each copy that you convey, and you may offer support or warranty protection for a fee.

5. Conveying Modified Source Versions.

You may convey a work based on the Program, or the modifications to produce it from the Program, in the form of source code under the terms of section 4, provided that you also meet all of these conditions:

- a) The work must carry prominent notices stating that you modified it, and giving a relevant date.
- b) The work must carry prominent notices stating that it is released under this License and any conditions added under section 7. This requirement modifies the requirement in section 4 to "keep intact all notices".
- c) You must license the entire work, as a whole, under this License to anyone who comes into possession of a copy. This License will therefore apply, along with any applicable section 7 additional terms, to the whole of the work, and all its parts, regardless of how they are packaged. This License gives no permission to license the work in any other way, but it does not invalidate such permission if you have separately received it.
- d) If the work has interactive user interfaces, each must display Appropriate Legal Notices; however, if the Program has interactive interfaces that do not display Appropriate Legal Notices, your work need not make them do so.

A compilation of a covered work with other separate and independent works, which are not by their nature extensions of the covered work, and which are not combined with it such as to form a larger program, in or on a volume of a storage or distribution medium, is called an "aggregate" if the compilation and its resulting copyright are not used to limit the access or legal rights of the compilation's users beyond what the individual works permit. Inclusion of a covered work in an aggregate does not cause this License to apply to the other parts of the aggregate.

6. Conveying Non-Source Forms.

You may convey a covered work in object code form under the terms of sections 4 and 5, provided that you also convey the machine-readable Corresponding Source under the terms of this License, in one of these ways:

- a) Convey the object code in, or embodied in, a physical product (including a physical distribution medium), accompanied by the Corresponding Source fixed on a durable physical medium customarily used for software interchange.

b) Convey the object code in, or embodied in, a physical product (including a physical distribution medium), accompanied by a written offer, valid for at least three years and valid for as long as you offer spare parts or customer support for that product model, to give anyone who possesses the object code either (1) a copy of the Corresponding Source for all the software in the product that is covered by this License, on a durable physical medium customarily used for software interchange, for a price no more than your reasonable cost of physically performing this conveying of source, or (2) access to copy the Corresponding Source from a network server at no charge.

c) Convey individual copies of the object code with a copy of the written offer to provide the Corresponding Source. This alternative is allowed only occasionally and noncommercially, and only if you received the object code with such an offer, in accord with subsection 6b.

d) Convey the object code by offering access from a designated place (gratis or for a charge), and offer equivalent access to the Corresponding Source in the same way through the same place at no further charge. You need not require recipients to copy the Corresponding Source along with the object code. If the place to copy the object code is a network server, the Corresponding Source may be on a different server (operated by you or a third party) that supports equivalent copying facilities, provided you maintain clear directions next to the object code saying where to find the Corresponding Source. Regardless of what server hosts the Corresponding Source, you remain obligated to ensure that it is available for as long as needed to satisfy these requirements.

e) Convey the object code using peer-to-peer transmission, provided you inform other peers where the object code and Corresponding Source of the work are being offered to the general public at no charge under subsection 6d.

A separable portion of the object code, whose source code is excluded from the Corresponding Source as a System Library, need not be included in conveying the object code work.

A "User Product" is either (1) a "consumer product", which means any tangible personal property which is normally used for personal, family, or household purposes, or (2) anything designed or sold for incorporation into a dwelling. In determining whether a product is a consumer product, doubtful cases shall be resolved in favor of coverage. For a particular product received by a particular user, "normally used" refers to a typical or common use of that class of product, regardless of the status of the particular user or of the way in which the particular user actually uses, or expects or is expected to use, the product. A product is a consumer product regardless of whether the product has substantial commercial, industrial or non-consumer uses, unless such uses represent the only significant mode of use of the product.

"Installation Information" for a User Product means any methods, procedures, authorization keys, or other information required to install and execute modified versions of a covered work in that User Product from a modified version of its Corresponding Source. The information must suffice to ensure that the continued functioning of the modified object code is in no case prevented or interfered with solely because modification has been made.

If you convey an object code work under this section in, or with, or specifically for use in, a User Product, and the conveying occurs as part of a transaction in which the right of possession and use of the User Product is transferred to the recipient in perpetuity or for a fixed term (regardless of how the transaction is characterized), the Corresponding Source conveyed under this section must be accompanied by the Installation Information. But this requirement does not apply if neither you nor any third party retains the ability to install modified object code on the User Product (for example, the work has been installed in ROM).

The requirement to provide Installation Information does not include a requirement to continue to provide support service, warranty, or updates for a work that has been modified or installed by the recipient, or for the User Product in which it has been modified or installed. Access to a network may be denied when the modification itself materially and adversely affects the operation of the network or violates the rules and protocols for communication across the network.

Corresponding Source conveyed, and Installation Information provided, in accord with this section must be in a format that is publicly documented (and with an implementation available to the public in source code form), and must require no special password or key for unpacking, reading or copying.

7. Additional Terms.

"Additional permissions" are terms that supplement the terms of this License by making exceptions from one or more of its conditions. Additional permissions that are applicable to the entire Program shall be treated as though they were included in this License, to the extent that they are valid under applicable law. If additional permissions apply only to part of the Program, that part may be used separately under those permissions, but the entire Program remains governed by this License without regard to the additional permissions.

When you convey a copy of a covered work, you may at your option remove any additional permissions from that copy, or from any part of it. (Additional permissions may be written to require their own removal in certain cases when you modify the work.) You may place additional permissions on material, added by you to a covered work, for which you have or can give appropriate copyright permission.

Notwithstanding any other provision of this License, for material you add to a covered work, you may (if authorized by the copyright holders of that material) supplement the terms of this License with terms:

- a) Disclaiming warranty or limiting liability differently from the terms of sections 15 and 16 of this License; or
- b) Requiring preservation of specified reasonable legal notices or author attributions in that material or in the Appropriate Legal Notices displayed by works containing it; or
- c) Prohibiting misrepresentation of the origin of that material, or requiring that modified versions of such material be marked in reasonable ways as different from the original version; or
- d) Limiting the use for publicity purposes of names of licensors or authors of the material; or
- e) Declining to grant rights under trademark law for use of some trade names, trademarks, or service marks; or
- f) Requiring indemnification of licensors and authors of that material by anyone who conveys the material (or modified versions of it) with contractual assumptions of liability to the recipient, for any liability that these contractual assumptions directly impose on those licensors and authors.

All other non-permissive additional terms are considered "further restrictions" within the meaning of section 10. If the Program as you received it, or any part of it, contains a notice stating that it is governed by this License along with a term that is a further restriction, you may remove that term. If a license document contains a further restriction but permits relicensing or conveying under this License, you may add to a covered work material governed by the terms of that license document, provided that the further restriction does not survive such relicensing or conveying.

If you add terms to a covered work in accord with this section, you must place, in the relevant source files, a statement of the additional terms that apply to those files, or a notice indicating where to find the applicable terms.

Additional terms, permissive or non-permissive, may be stated in the form of a separately written license, or stated as exceptions; the above requirements apply either way.

8. Termination.

You may not propagate or modify a covered work except as expressly provided under this License. Any attempt otherwise to propagate or modify it is void, and will automatically terminate your rights under this License (including any patent licenses granted under the third paragraph of section 11).

However, if you cease all violation of this License, then your license from a particular copyright holder is reinstated (a) provisionally, unless and until the copyright holder explicitly and finally terminates your license, and (b) permanently, if the copyright holder fails to notify you of the violation by some reasonable means prior to 60 days after the cessation.

Moreover, your license from a particular copyright holder is reinstated permanently if the copyright holder notifies you of the violation by some reasonable means, this is the first time you have received notice of violation of this License (for any work) from that copyright holder, and you cure the violation prior to 30 days after your receipt of the notice.

Termination of your rights under this section does not terminate the licenses of parties who have received copies or rights from you under this License. If your rights have been terminated and not permanently reinstated, you do not qualify to receive new licenses for the same material under section 10.

9. Acceptance Not Required for Having Copies.

You are not required to accept this License in order to receive or run a copy of the Program. Ancillary propagation of a covered work occurring solely as a consequence of using peer-to-peer transmission to receive a copy likewise does not require acceptance. However, nothing other than this License grants you permission to propagate or modify any covered work. These actions infringe copyright if you do not accept this License. Therefore, by modifying or propagating a covered work, you indicate your acceptance of this License to do so.

10. Automatic Licensing of Downstream Recipients.

Each time you convey a covered work, the recipient automatically receives a license from the original licensors, to run, modify and propagate that work, subject to this License. You are not responsible for enforcing compliance by third parties with this License.

An "entity transaction" is a transaction transferring control of an organization, or substantially all assets of one, or subdividing an organization, or merging organizations. If propagation of a covered work results from an entity transaction, each party to that transaction who receives a copy of the work also receives whatever licenses to the work the party's predecessor in interest had or could give under the previous paragraph, plus a right to possession of the Corresponding Source of the work from the predecessor in interest, if the predecessor has it or can get it with reasonable efforts.

You may not impose any further restrictions on the exercise of the rights granted or affirmed under this License. For example, you may not impose a license fee, royalty, or other charge for exercise of rights granted under this License, and you may not initiate litigation (including a cross-claim or counterclaim in a lawsuit) alleging that any patent claim is infringed by making, using, selling, offering for sale, or importing the Program or any portion of it.

11. Patents.

A "contributor" is a copyright holder who authorizes use under this License of the Program or a work on which the Program is based. The work thus licensed is called the contributor's "contributor version".

A contributor's "essential patent claims" are all patent claims owned or controlled by the contributor, whether already acquired or hereafter acquired, that would be infringed by some manner, permitted by this License, of making, using, or selling its contributor version, but do not include claims that would be infringed only as a consequence of further modification of the contributor version. For purposes of this definition, "control" includes the right to grant patent sublicenses in a manner consistent with the requirements of this License.

Each contributor grants you a non-exclusive, worldwide, royalty-free patent license under the contributor's essential patent claims, to make, use, sell, offer for sale, import and otherwise run, modify and propagate the contents of its contributor version.

In the following three paragraphs, a "patent license" is any express agreement or commitment, however denominated, not to enforce a patent (such as an express permission to practice a patent or covenant not to sue for patent infringement). To "grant" such a patent license to a party means to make such an agreement or commitment not to enforce a patent against the party.

If you convey a covered work, knowingly relying on a patent license, and the Corresponding Source of the work is not available for anyone to copy, free of charge and under the terms of this License, through a publicly available network server or other readily accessible means, then you must either (1) cause the Corresponding Source to be so available, or (2) arrange to deprive yourself of the benefit of the patent license for this particular work, or (3) arrange, in a manner consistent with the requirements of this License, to extend the patent license to downstream recipients. "Knowingly relying" means you have actual knowledge that, but for the patent license, your conveying the covered work in a country, or your recipient's use of the covered work in a country, would infringe one or more identifiable patents in that country that you have reason to believe are valid.

If, pursuant to or in connection with a single transaction or arrangement, you convey, or propagate by procuring conveyance of, a covered work, and grant a patent license to some of the parties receiving the covered work authorizing them to use, propagate, modify or convey a specific copy of the covered work, then the patent license you grant is automatically extended to all recipients of the covered work and works based on it.

A patent license is "discriminatory" if it does not include within the scope of its coverage, prohibits the exercise of, or is conditioned on the non-exercise of one or more of the rights that are specifically granted under this License. You may not convey a covered work if you are a party to an arrangement with a third party that is in the business of distributing software, under which you make payment to the third party based on the extent of your activity of conveying the work, and under which the third party grants, to any of the parties who would receive the covered work from you, a discriminatory patent license (a) in connection with copies of the covered work conveyed by you (or copies made from those copies), or (b) primarily for and in connection with specific products or compilations that contain the covered work, unless you entered into that arrangement, or that patent license was granted, prior to 28 March 2007.

Nothing in this License shall be construed as excluding or limiting any implied license or other defenses to infringement that may otherwise be available to you under applicable patent law.

12. No Surrender of Others' Freedom.

If conditions are imposed on you (whether by court order, agreement or otherwise) that contradict the conditions of this License, they do not excuse you from the conditions of this License. If you cannot convey a covered work so as to satisfy simultaneously your obligations under this License and any other pertinent obligations, then as a consequence you may not convey it at all. For example, if you agree to terms that obligate you to collect a royalty for further conveying from those to whom you convey the Program, the only way you could satisfy both those terms and this License would be to refrain entirely from conveying the Program.

13. Use with the GNU Affero General Public License.

Notwithstanding any other provision of this License, you have permission to link or combine any covered work with a work licensed under version 3 of the GNU Affero General Public License into a single combined work, and to convey the resulting work. The terms of this License will continue to apply to the part which is the covered work, but the special requirements of the GNU Affero General Public License, section 13, concerning interaction through a network will apply to the combination as such.

14. Revised Versions of this License.

The Free Software Foundation may publish revised and/or new versions of the GNU General Public License from time to time. Such new versions will be similar in spirit to the present version, but may differ in detail to address new problems or concerns.

Each version is given a distinguishing version number. If the Program specifies that a certain numbered version of the GNU General Public License "or any later version" applies to it, you have the option of following the terms and conditions either of that numbered version or of any later version published by the Free Software Foundation. If the Program does not specify a version number of the GNU General Public License, you may choose any version ever published by the Free Software Foundation.

If the Program specifies that a proxy can decide which future versions of the GNU General Public License can be used, that proxy's public statement of acceptance of a version permanently authorizes you to choose that version for the Program.

Later license versions may give you additional or different permissions. However, no additional obligations are imposed on any author or copyright holder as a result of your choosing to follow a later version.

15. Disclaimer of Warranty.

THERE IS NO WARRANTY FOR THE PROGRAM, TO THE EXTENT PERMITTED BY APPLICABLE LAW. EXCEPT WHEN OTHERWISE STATED IN WRITING THE COPYRIGHT HOLDERS AND/OR OTHER PARTIES PROVIDE THE PROGRAM "AS IS" WITHOUT WARRANTY OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE ENTIRE RISK AS TO THE QUALITY AND PERFORMANCE OF THE PROGRAM IS WITH YOU. SHOULD THE PROGRAM PROVE DEFECTIVE, YOU ASSUME THE COST OF ALL NECESSARY SERVICING, REPAIR OR CORRECTION.

16. Limitation of Liability.

IN NO EVENT UNLESS REQUIRED BY APPLICABLE LAW OR AGREED TO IN WRITING WILL ANY COPYRIGHT HOLDER, OR ANY OTHER PARTY WHO MODIFIES AND/OR CONVEYS THE PROGRAM AS PERMITTED ABOVE, BE LIABLE TO YOU FOR DAMAGES, INCLUDING ANY GENERAL, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF THE USE OR INABILITY TO USE THE PROGRAM (INCLUDING BUT NOT LIMITED TO LOSS OF DATA OR DATA BEING RENDERED INACCURATE OR LOSSES SUSTAINED BY YOU OR THIRD PARTIES OR A FAILURE OF THE PROGRAM TO OPERATE WITH ANY OTHER PROGRAMS), EVEN IF SUCH HOLDER OR OTHER PARTY HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

17. Interpretation of Sections 15 and 16.

If the disclaimer of warranty and limitation of liability provided above cannot be given local legal effect according to their terms, reviewing courts shall apply local law that most closely approximates an absolute waiver of all civil liability in connection with the Program, unless a warranty or assumption of liability accompanies a copy of the Program in return for a fee.

END OF TERMS AND CONDITIONS

How to Apply These Terms to Your New Programs

If you develop a new program, and you want it to be of the greatest possible use to the public, the best way to achieve this is to make it free software which everyone can redistribute and change under these terms.

To do so, attach the following notices to the program. It is safest to attach them to the start of each source file to most effectively state the exclusion of warranty; and each file should have at least the "copyright" line and a pointer to where the full notice is found.

<one line to give the program's name and a brief idea of what it does.>

Copyright (C) <year> <name of author>

This program is free software: you can redistribute it and/or modify it under the terms of the GNU General Public License as published by the Free Software Foundation, either version 3 of the License, or (at your option) any later version.

This program is distributed in the hope that it will be useful, but WITHOUT ANY WARRANTY; without even the implied warranty of MERCHANTABILITY or FITNESS FOR A PARTICULAR PURPOSE. See the GNU General Public License for more details.

You should have received a copy of the GNU General Public License along with this program. If not, see <<http://www.gnu.org/licenses/>>.

Also add information on how to contact you by electronic and paper mail.

If the program does terminal interaction, make it output a short notice like this when it starts in an interactive mode:

<program> Copyright (C) <year> <name of author>

This program comes with ABSOLUTELY NO WARRANTY; for details type `show w'.

This is free software, and you are welcome to redistribute it under certain conditions; type `show c' for details.

The hypothetical commands `show w' and `show c' should show the appropriate parts of the General Public License. Of course, your program's commands might be different; for a GUI interface, you would use an "about box".

You should also get your employer (if you work as a programmer) or school, if any, to sign a "copyright disclaimer" for the program, if necessary.

For more information on this, and how to apply and follow the GNU GPL, see <<http://www.gnu.org/licenses/>>.

The GNU General Public License does not permit incorporating your program into proprietary programs. If your program is a subroutine library, you may consider it more useful to permit linking proprietary applications with the library. If this is what you want to do, use the GNU Lesser General Public License instead of this License. But first, please read <<http://www.gnu.org/philosophy/why-not-lgpl.html>>.

VII. The LibPng License

This copy of the libpng notices is provided for your convenience. In case of any discrepancy between this copy and the notices in the file png.h that is included in the libpng distribution, the latter shall prevail.

COPYRIGHT NOTICE, DISCLAIMER, and LICENSE:

If you modify libpng you may insert additional notices immediately following this sentence.

libpng versions 1.2.6, August 15, 2004, through 1.2.36, May 7, 2009, are Copyright (c) 2004, 2006-2009 Glenn Randers-Pehrson, and are distributed according to the same disclaimer and license as libpng-1.2.5 with the following individual added to the list of Contributing Authors

Cosmin Truta

libpng versions 1.0.7, July 1, 2000, through 1.2.5 - October 3, 2002, are Copyright (c) 2000-2002 Glenn Randers-Pehrson, and are distributed according to the same disclaimer and license as libpng-1.0.6 with the following individuals added to the list of Contributing Authors

Simon-Pierre Cadieux

Eric S. Raymond

Gilles Vollant

and with the following additions to the disclaimer:

There is no warranty against interference with your enjoyment of the library or against infringement. There is no warranty that our efforts or the library will fulfill any of your particular purposes or needs. This library is provided with all faults, and the entire risk of satisfactory quality, performance, accuracy, and effort is with the user.

libpng versions 0.97, January 1998, through 1.0.6, March 20, 2000, are Copyright (c) 1998, 1999 Glenn Randers-Pehrson, and are distributed according to the same disclaimer and license as libpng-0.96, with the following individuals added to the list of Contributing Authors:

Tom Lane

Glenn Randers-Pehrson

Willem van Schaik

libpng versions 0.89, June 1996, through 0.96, May 1997, are Copyright (c) 1996, 1997 Andreas Dilger Distributed according to the same disclaimer and license as libpng-0.88, with the following individuals added to the list of Contributing Authors:

John Bowler
Kevin Bracey
Sam Bushell
Magnus Holmgren
Greg Roelofs
Tom Tanner

libpng versions 0.5, May 1995, through 0.88, January 1996, are Copyright (c) 1995, 1996 Guy Eric Schalnat, Group 42, Inc.

For the purposes of this copyright and license, "Contributing Authors" is defined as the following set of individuals:

Andreas Dilger
Dave Martindale
Guy Eric Schalnat
Paul Schmidt
Tim Wegner

The PNG Reference Library is supplied "AS IS". The Contributing Authors and Group 42, Inc. disclaim all warranties, expressed or implied, including, without limitation, the warranties of merchantability and of fitness for any purpose. The Contributing Authors and Group 42, Inc. assume no liability for direct, indirect, incidental, special, exemplary, or consequential damages, which may result from the use of the PNG Reference Library, even if advised of the possibility of such damage.

Permission is hereby granted to use, copy, modify, and distribute this source code, or portions hereof, for any purpose, without fee, subject to the following restrictions:

1. The origin of this source code must not be misrepresented.
2. Altered versions must be plainly marked as such and must not be misrepresented as being the original source.
3. This Copyright notice may not be removed or altered from any source or altered source distribution.

The Contributing Authors and Group 42, Inc. specifically permit, without fee, and encourage the use of this source code as a component to supporting the PNG file format in commercial products. If you use this source code in a product, acknowledgment is not required but would be appreciated.

A "png_get_copyright" function is available, for convenient use in "about" boxes and the like:

```
printf("%s",png_get_copyright(NULL));
```

Also, the PNG logo (in PNG format, of course) is supplied in the files "pngbar.png" and "pngbar.jpg" (88x31) and "pngnow.png" (98x31).

Libpng is OSI Certified Open Source Software. OSI Certified Open Source is a certification mark of the Open Source Initiative.

Glenn Randers-Pehrson
glennrp at users.sourceforge.net
May 7, 2009

VIII. The zlib License

Copyright (c) <"year"> <"copyright holders">

This software is provided 'as-is', without any express or implied warranty. In no event will the authors be held liable for any damages arising from the use of this software.

Permission is granted to anyone to use this software for any purpose, including commercial applications, and to alter it and redistribute it freely, subject to the following restrictions:

1. The origin of this software must not be misrepresented; you must not claim that you wrote the original software. If you use this software in a product, an acknowledgment in the product documentation would be appreciated but is not required.
2. Altered source versions must be plainly marked as such, and must not be misrepresented as being the original software.
3. This notice may not be removed or altered from any source distribution.

